

## Compound Digest Solution

<b>Catalog No.</b>	AR0022
<b>Size</b>	50ml
<b>Type</b>	Ready to use
<b>Storage</b>	Store at 4°C for one year.

### Product Description

Compound Digest Solution was manufactured by Boster, it is used for enzymic digestion in immunohistochemistry (IHC). It contains enzymes that catalyze the splitting of proteins into smaller peptide fractions and amino acids by a process known as proteolysis.

### Main Contents

Multiple proteinases

### Application

1. Add Compound Digest Solution directly on to the tissue or cell, and cover the tissue or cell completely,
2. Incubate for 5-10 minutes in 37°C incubator,
3. Take out the tissue or cell and wash it with PBS.