



17 α 20 β -diOH-Progesteron-3-CMO-BSA polyclonal antibody

Product No. FKA-332

The Cross reactivities of 17 α 20 β -diOH-Progesteron-CMO-BSA (FKA-332) by RI Assay

Antigen	17 α 20 β -diOH-Progesteron-CMO-BSA
Antibody[FKA-332]	lyophilized Serum 5 μ l / Vial
Dilution of Serum	2000test for RI Assay 250 μ l / tube use
Binding of Total Testosterone dose of B0/2	
Strage/Handling	Maintain at 4 $^{\circ}$ C-8 $^{\circ}$ C (-20 $^{\circ}$ C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

Steroid	Cross Reaction
17 α 20 β - Progesteron	100.00 %
17 α 20 α - Progesteron	3.40 %
17 α -OH- Progesteron	1.00 %
Cortisol	0.01 %

For research use only; not for use as a diagnostic.

