# Homo sapiens DCTD cDNA Clone

Catalog Number: HG14664-G



Gene : dCMP deaminase

Official Syr	nbol :	DCTD
--------------	--------	------

- Synonym : DCTD
- Source : Homo sapiens
- cDNA Size: 537
- **RefSeq:** BC001286

## Description

Lot : Please refer to the label on the tube

## **Sequence Description :**

Identical with the Gene Bank Ref. ID sequence.

Vector:

pGEM-T

Shipping carrier :

Each tube contains approximately 10 µg of lyophilized plasmid.

## Storage :

The lyophilized plasmid can be stored at ambient temperature for three months.

## **Quality control:**

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

## Sequencing primer list :

M13-47 : 5' GCCAGGGTTTTCCCAGTCACGAC 3'

RV-M: 5' GAGCGGATAACAATTTCACACAGG 3'

Other M13 primers can also be used as sequencing primers.



## **Plasmid Resuspension protocol**

1.Centrifuge the tube for 5~10 min at 4,000 rpm.

2.Carefully open the tube and add 100  $\mu$ l of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature. 4.Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 4000 rpm. 5.Store the plasmid at -20 °C.

#### The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

## E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. DH5 α, TOP10, JM109.


# Homo sapiens DCTD cDNA Clone

Catalog Number: HG14664-G

## Vector Information

The pGEM-T vector is a high-efficiency TA cloning vector which contains multiple cloning sites as shown below. The pGEM-T vector is 3.0kb in size and contains the amplicin resistance gene for selection. The coding sequence was inserted by TA cloning.

#### Physical Map of pGEM-T :



 Please refer to http://www.sinobiological.com/Vector-pGEM-T-a-1636.html for the vector sequence.

