# Homo sapiens FAM19A3 cDNA Clone

Catalog Number: HG15860-G

## **General Information**

 Gene :
 family with sequence similarity 19 (chemokine (C-C motif)-like), member A3

 Official Symbol :
 FAM19A3

Synonym : TAFA3, TAFA-3

- Source : Homo sapiens
- **cDNA Size:** 510
- RefSeq : BC112268

# Description

Lot : Please refer to the label on the tube

**Sequence Description :** 

Identical with the Gene Bank Ref. ID sequence.

Vector :

pGEM-T

Shipping carrier :

Each tube contains approximately 10 µg of lyophilized plasmid.

#### Storage :

The lyophilized plasmid can be stored at ambient temperature for three months.

#### Quality control :

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

#### Sequencing primer list :

RV-M : 5' GAGCGGATAACAATTTCACACAGG 3'

Other M13 primers can also be used as sequencing primers.



# **Plasmid Resuspension protocol**

1.Centrifuge the tube for 5~10 min at 4,000 rpm.

2.Carefully open the tube and add 100  $\mu l$  of sterile water to dissolve the DNA.

3.Close the tube and incubate for 10 minutes at room temperature. 4.Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 4000 rpm. 5.Store the plasmid at -20  $^{\circ}$ C.

### The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

# *E.coli* strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. DH5  $\alpha,$  TOP10, JM109.


# Homo sapiens FAM19A3 cDNA Clone

Catalog Number: HG15860-G

## Vector Information

The pGEM-T vector is a high-efficiency TA cloning vector which contains multiple cloning sites as shown below. The pGEM-T vector is 3.0kb in size and contains the amplicin resistance gene for selection. The coding sequence was inserted by TA cloning.

### Physical Map of pGEM-T :



 Please refer to http://www.sinobiological.com/Vector-pGEM-T-a-1636.html for the vector sequence.

