
Product Manual


CytoSelect™ Leukocyte- Endothelium Adhesion Assay

Catalog Number

CBA- 210

100 assays

FOR RESEARCH USE ONLY
Not for use in diagnostic procedures



CELL BIOLABS, INC.

Creating Solutions for Life Science Research

Introduction

Leukocyte extravasation into perivascular tissue plays a key role in inflammatory diseases. This recruitment requires leukocyte interaction with vascular endothelium and consists of multiple, consecutive processes including the capture of circulating leukocytes, subsequent leukocyte rolling, arrest, firm adhesion and transmigration (Figure 1). This multistep paradigm is realized by sequential activation-dependent interactions between endothelial cell adhesion molecules and their specific ligands on leukocytes. The first step of transient adhesion and rolling is known to be mediated by an interaction of leukocyte or endothelial cell selectins and their oligosaccharide-bearing ligands. Arrest and firm adhesion of leukocytes to endothelium is dependent on the activation of $\beta 2$ integrins like Mac-1 or LFA-1 on the leukocyte cell surface, followed by interaction with endothelial cell proteins belonging to the Ig superfamily such as ICAM-1.

Cell Biolabs' CytoSelect™ Leukocyte-endothelium Adhesion Assay provides a robust system for the quantitative determination of leukocyte-endothelium interactions. The kit contains sufficient reagents for the evaluation of 100 assays in a 96-well plate.





Figure 1. The Leukocyte-Endothelium Adhesion Cascade.

Assay Principle



Related Products

1. CBA-052: CytoSelect™ 48-Well Cell Adhesion Assay (Collagen I, Colorimetric)
2. CBA-053: CytoSelect™ 48-Well Cell Adhesion Assay (Collagen I, Fluorometric)
3. CBA-070: CytoSelect™ 48-Well Cell Adhesion Assay (ECM Array, Colorimetric)
4. CBA-071: CytoSelect™ 48-Well Cell Adhesion Assay (ECM Array, Fluorometric)
5. CBA-200: Endothelial Tube Formation Assay (In Vitro Angiogenesis)
6. CBA-211: CytoSelect™ Leukocyte-Epithelium Adhesion Assay
7. CBA-215: CytoSelect™ Tumor-Endothelium Adhesion Assay
8. CBA-212: CytoSelect™ Leukocyte Transmigration Assay
9. CBA-216: CytoSelect™ Tumor Transendothelial Migration Assay
10. CBA-320: CytoSelect™ 96-Well Hematopoietic Colony Forming Cell Assay

Kit Components

1. 500X LeukoTracker™ Solution (Part No. 12101): One 100 μL tube
2. Gelatin Solution (Part No. 12102): One 12 mL bottle of sterile 0.1 % Gelatin in 1X PBS
3. 4X Lysis Buffer (Part No. 10404): One 10 mL bottle

4. 10X Wash Buffer (Part No. 12104): One 20 mL bottle
5. TNF α (Part No. 12105): One 100 μ L tube of 10 μ g/mL TNF α in sterile 1X PBS/0.1%BSA

Materials Not Supplied

1. Endothelial cells and cell culture medium
2. Leukocytes
3. 96-well or 48-well tissue culture plate
4. Serum free medium, such as DMEM containing 0.5% BSA, 2 mM CaCl₂ and 2 mM MgCl₂
5. Sterile 1X PBS
6. Cell culture incubator (37°C, 5% CO₂ atmosphere)
7. Light microscope
8. 96-well plate suitable for a fluorescence plate reader
9. Fluorescence plate reader

Storage

LeukoTracker™ Solution and TNF α should be removed from the kit and stored at -20°C immediately. Store all other components at 4°C.

Preparation of Reagents

- 1X Wash Buffer: Prepare a 1X Wash Buffer by diluting the provided 10X stock 1:10 in deionized water. Store the diluted solution at room temperature.
- 1X Lysis Buffer: Prepare a 1X Lysis Buffer by diluting the provided 4X stock 1:4 in deionized water. Store the diluted solution at room temperature.

Gelatin Coating

1. Under sterile conditions, add 200 μ L of the Gelatin Solution to each well of a 48-well tissue culture treated plate, or 100 μ L of the Gelatin Solution to each well of a 96-well tissue culture treated plate.
2. Incubate for 60 min at 37°C in a cell culture incubator.
3. Wash twice with sterile 1X PBS. Aspirate the final wash before use.

Assay Protocol

1. Add 50,000-100,000 endothelial cells/well to the Gelatin-coated 48-well or 96-well plate.
2. Culture cells for 48-72 until the endothelial cells form a monolayer.
3. Treat endothelial cell monolayer or leukocyte with desired activator or inhibitor for 6-12 hrs.

4. Harvest leukocytes and prepare a cell suspension at 1.0×10^6 cells/ml in serum free media. Add LeukoTracker to a final concentration of 1X (for example, add 2 μ L of 500X LeukoTracker™ solution to 1.0 mL of leukocyte cell suspension).
5. Incubate for 60 min at 37°C in a cell culture incubator. Spin down cells at 1000 rpm for 2 minutes, aspirate the medium and wash cell pellet with serum free media. Repeat the wash twice. Resuspend the cell pellet at $0.25 - 1.0 \times 10^6$ cells/ml in serum free media.
6. Aspirate endothelial culture media and wash once with serum free media. Add 200 μ L of the cell suspension to each well already containing the endothelial monolayer.
7. Incubate for 30-90 min in a cell culture incubator.
8. **Carefully** discard or aspirate the media from each well (**Note: Do not allow wells to dry**). Gently wash each well 3 times with 250 μ L 1X Wash Buffer.
9. (Optional) Count the adherent leukocytes under an inverted fluorescence microscope; average at least three separate fields per well.
10. Aspirate the final wash and add 150 μ L of 1X Lysis Buffer to each well containing cells. Incubate 5 minutes at room temperature with shaking.
11. Transfer 100 μ L of the mixture to a 96-well plate suitable for fluorescence measurement. Read fluorescence with a fluorescence plate reader at 480 nm/520 nm.

Example of Results

The following figures demonstrate typical with Cell Biolabs CytoSelect™ Leukocyte-endothelium Adhesion Assay Kit. Fluorescence measurement was performed on SpectraMax Gemini XS Fluorometer (Molecular Devices) with a 485/538 nm filter set and 530 nm cutoff. One should use the data below for reference only. This data should not be used to interpret actual results.




Figure 2. Quantitation of Human Monocytic THP-1. LeukoTracker™ labeled THP-1 cells were titrated in 1X PBS, then subsequently lysed with 2X Lysis Buffer (75 μ L of cell suspension was mixed with 75 μ L of 2X Lysis Buffer). Fluorescence was quantified as described in Assay Protocol.




Figure 3. Human Monocytic THP-1 Adhesion to HUVEC Monolayer. HUVEC monolayer in 48-well plate was treated with 1 μ M PMA for 12 hrs. LeukoTracker™ labeled THP-1 cells were allowed to attach to HUVEC monolayer for 1 hr. Adherent cells were lysed and quantified by as described in the Assay Protocol.





Figure 4. Cytokine Activation of Human Monocytic THP-1 Adhesion to HUVEC. HUVEC monolayer (left) or THP-1 cells (right) were treated with 50 ng/mL TNF α or 1 μ M PMA for 12 hrs. LeukoTracker™ labeled THP-1 cells (50,000 cells/well) were allowed to attach to HUVEC monolayer for 1 hr. Adherent cells were lysed and quantified by as described in the Assay Protocol.

References

1. Carlos T. M and Harlan J. M. (1994) *Blood* **84**, 2068-101.
2. Meldrum, D. R. (1998) *Am. J. Physiol.* **274**, R577-R595.

Recent Product Citations

1. Wojkowska, D.J. et al. (2017). Interleukin 17A Promotes Lymphocytes Adhesion and Induces CCL2 and CXCL1 Release from Brain Endothelial Cells. *Int. J. Mol. Sci.* doi:10.3390/ijms18051000.
2. Wojkowska, D. et al. (2017). Th17 Cells Interactions with the Brain Endothelium in Vitro. *Preprints*. 2017030094. doi: 10.20944/preprints201703.0094.v1.
3. He, M. et al. (2016). miR-483 targeting of CTGF suppresses endothelial-to-mesenchymal transition: therapeutic implications in Kawasaki Disease. *Circ. Res.* **120**(2):354-365. doi:10.1161/CIRCRESAHA.116.310233.
4. Rossignoli A, et al. (2017). Poliovirus Receptor-Related 2: A Cholesterol-Responsive Gene Affecting Atherosclerosis Development by Modulating Leukocyte Migration. *Arterioscler Thromb Vasc Biol.* **37**(3):534-542. doi: 10.1161/ATVBAHA.116.308715.
5. Cao, Q. et al. (2016). Myeloid deletion of α 1AMPK exacerbates atherosclerosis in LDL receptor knockout (LDLRKO) mice. *Diabetes.* **65**:1565-1576.
6. Huang, M. et al. (2015). Niclosamide inhibits the inflammatory and angiogenic activation of human umbilical vein endothelial cells. *Inflamm Res.* **64**:1023-1032.
7. Shah, D. et al. (2015). C1q deficiency promotes pulmonary vascular inflammation and enhances the susceptibility of the lung endothelium to injury. *J Biol Chem.* doi:10.1074/jbc.M115.690784.
8. Campos-Estrada, C. et al. (2015). Simvastatin and benznidazole-mediated prevention of Trypanosoma cruzi-induced endothelial activation: role of 15-epi-lipoxin A4 in the action of simvastatin. **9**:e0003770.
9. Ibrahim, A. S. et al. (2015). A lipidomic screen of hyperglycemia-treated HRECs links 12/15-Lipoxygenase to microvascular dysfunction during diabetic retinopathy via NADPH oxidase. *J Lipid Res.* **56**:599-611.
10. Cao, Q. et al. (2014). Inhibiting DNA Methylation by 5-Aza-2'-deoxycytidine ameliorates atherosclerosis through suppressing macrophage inflammation. *Endocrinology.* **155**:4925-4938.
11. Cary, L. H. et al. (2014). Interactions between endothelial cells and T cells modulate responses to mixed neutron/gamma radiation. *Radiat Res.* **181**:592-604.
12. Burrows, G. G. et al. (2014). Human multipotent adult progenitor cells transcriptionally regulate fucosyltransferase VII. *Cytotherapy.* **16**:566-575.
13. Kapitsinou, P. P. et al. (2014). Endothelial HIF-2 mediates protection and recovery from ischemic kidney injury. *J Clin Invest.* **124**:2396.
14. Wu, X. Y. et al. (2014). Regulation of microRNA-155 in endothelial inflammation by targeting nuclear factor (NF)- κ B p65. *J Cell Biochem.* **115**:1928-1936.
15. Ghoshal, P. et al. (2014). Glycosylation inhibitors efficiently inhibit P-selectin-mediated cell adhesion to endothelial cells. *PLoS One.* **9**:e99363.
16. Sharma-Walia, N. et al. (2014). The Kaposi's Sarcoma-associated Herpesvirus (KSHV)-induced 5-lipoxygenase-leukotriene B4 cascade plays key roles in KSHV latency, monocyte recruitment, and lipogenesis. *J. Virol.* **88**:2131-2156.
17. Avan, A. et al. (2013). Norepinephrine increases NADPH oxidase-derived superoxide in human peripheral blood mononuclear cells via α -adrenergic receptors. *Cancer Res.* **73**:6745-6756

18. Liu, G. et al. (2012). ICAM-1–activated Src and eNOS signaling increase endothelial cell surface PECAM-1 adhesivity and neutrophil transmigration. *Blood*. **120**:1942-1952.
19. Fang, Y. et al. (2012). Site-Specific microRNA-92a regulation of Kruppel-like factors 4 and 2 in atherosusceptible endothelium. *Arterioscler Thromb Vasc Biol*. **32**:979-987.
20. Curatola, A.M. et al. (2012). Dehydroepiandrosterone (DHEA) inhibition of monocyte binding by vascular endothelium is associated with sialylation of neural cell adhesion molecule. *Reproductive Sciences* **19**:86-91.
21. Wang, Y.L. et al. (2011). Innate immune function of the adherens junction protein p120-Catenin in endothelial response to endotoxin. *J. Immunol*. **186**-3180-3187.
22. Xue, J. et al. (2009). NFkB regulates thrombin-induced ICAM-1 gene expression in cooperation with NFAT by binding to the intronic NFkB Site in the ICAM-1 gene. *Physiol. Genomics* 10.1152/physiolgenomics.00012.2009.

Warranty

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS' sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

Contact Information

Cell Biolabs, Inc.
7758 Arjons Drive
San Diego, CA 92126
Worldwide: +1 858-271-6500
USA Toll-Free: 1-888-CBL-0505
E-mail: tech@cellbiolabs.com
www.cellbiolabs.com

©2004-2017: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.