



Anti Estradiol-3-carboxymethylether-BSA polyclonal antibody

Product No. FKA-236E

The Cross reactivities of anti Estradiol-3-carboxymethylether-BSA (FKA-236E)
by ELISA with Estradiol-3-carboxymethylethe-HRP (FKA-235)

Antigen	Estradiol-3-carboxymethylether-BSA		
Antibody(FKA-236E)	Purified by Protein A	lyophilized	2μl / Vial
Dilution of IgG	2000test for ELISA	50μl/well use	
Label (FKA-235)	Estradiol-3-carboxymethylethe-HRP		1ml / Vial
	1,000test for ELISA	50μl/well use	
O.D. (492/600nm) of B ₀	about 0.9		
Estradiol dose of B ₀ /2	4 ng/ ml		
Strage/Handling	Maintain at 4°C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.		

Steroid	Cross Reaction
Estradiol(E2)	100 %
Estrone(E1)	0.80 %
Estriol(E3)	0.50 %
Estrone-3-Sulfate	26.80 %
Estradiol-3-glucuronide	56.30 %
Estrone-3-glucuronide	1.20 %
Testosterone	0.05 %
Pregnenolone	0.00 %
Progesterone	0.00 %
Cortisol	0.00 %
Cortisone	0.00 %
4-Androstenedione	0.00 %
DHA	0.00 %
DHT	0.00 %

For research use only; not for use as a diagnostic.

