



**MONOCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. 2ZHIGA2**

# Anti Human IgA2

<b>Product type</b>	Primary antibodies
<b>Host</b>	Mouse
<b>Source</b>	Ascites Hybridoma cell line (Hyb-7311A) was grown in ascitic fluid BALB/c Mouse
<b>Form</b>	Liquid Ammonium sulfate precipitation Physiological salt solution containing 0.1% Sodium Azide
<b>Volume</b>	0.5 mL
<b>Concentration</b>	0.5 mg/mL
<b>Specificity</b>	This antibody specifically reacts with Human IgA2
<b>Antigen</b>	IgA1 loss syndrome patient origin IgA2 antigen
<b>Clone</b>	Hyb-7311A The splenic lymphocytes from BALB/c mouse immunized with IgA1 loss syndrome patient origin IgA2 were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-7311A) were screened.
<b>Isotype</b>	IgG <sub>1</sub>

**Application notes** ELISA

**Recommended use**

**Recommended dilutions**

Optimal dilutions/concentrations should be determined by the end user.

**Staining Pattern**

**Cross reactivity**

**Storage**

Maintain at 2-8 °C (-20 °C for long term storage) in undiluted aliquots.  
Avoid repeated freeze/thaw cycles.

**References**

1) Hikata M, et al: Hepatology 6(4): 652-657, 1986.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619



## 抗ヒト IgA2

<b>Product type</b>	一次抗体
<b>Host</b>	マウス
<b>Source</b>	BALB/c マウス腹水 ハイブリドーマ (Hyb-7311A) を BALB/c マウスに接種して得られた腹水
<b>Form</b>	液状 硫酸アンモニウム塩析精製 0.1%アジ化ナトリウムを含む生理的食塩水溶液
<b>Volume</b>	0.5 mL
<b>Concentration</b>	0.5 mg/mL
<b>Specificity</b>	ヒト IgA2 に特異的に反応
<b>Antigen</b>	IgA1 欠損症患者血清由来 IgA2 抗原
<b>Clone</b>	Hyb-7311A IgA1 欠損症患者血清由来 IgA2 を免疫した BALB/c マウスの脾臓細胞とマウスミエローマ細胞 (NS-1) を融合して得られたハイブリドーマ
<b>Isotype</b>	IgG <sub>1</sub>

**Application notes**

ELISA

**Recommended use**

**Recommended dilutions**

Optimal dilutions/concentrations should be determined by the end user.

**Staining Pattern**

**Cross reactivity**

**Storage**

2~8°C                      3ヶ月間

-20°C                        3年間

(凍結融解の繰り返しは避けてください)

**References**

1) Hikata M, et al: Hepatology 6(4): 652-657, 1986.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619