



MONOCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. 2ZHIGE1

Anti Human IgE

Product type	Primary antibodies
Host	Mouse
Source	Ascites Hybridoma cell line (Hyb-IgE-4) was grown in ascitic fluid BALB/c Mouse
Form	Liquid Ammonium sulfate precipitation Physiological salt solution containing 0.1% Sodium Azide
Volume	0.5 mL
Concentration	0.5 mg/mL
Specificity	This antibody specifically reacts with Human IgE
Antigen	Purified Human IgE antigen
Clone	Hyb-IgE-4 The splenic lymphocytes from BALB/c mouse immunized with purified Human IgE were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-IgE-4) were screened.
Isotype	IgG ₁

Application notes ELISA
Recommended use

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Storage Maintain at 2-8 °C (-20 °C for long term storage) in undiluted aliquots.
Avoid repeated freeze/thaw cycles.

References



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618 FAX : +81-3-5632-9619



モノクローナル抗体

For research use only. Not for clinical diagnosis.

Catalog No. 2ZHIGE1

抗ヒト IgE

Product type	一次抗体
Host	マウス
Source	BALB/c マウス腹水 ハイブリドーマ (Hyb-IgE-4) を BALB/c マウスに接種して得られた腹水
Form	液状 硫酸アンモニウム塩析精製 0.1%アジ化ナトリウムを含む生理的食塩水溶液
Volume	0.5 mL
Concentration	0.5 mg/mL
Specificity	ヒト IgE に特異的に反応
Antigen	精製ヒト IgE 抗原
Clone	Hyb-IgE-4 精製ヒト IgE を免疫した BALB/c マウスの脾臓細胞とマウスミエローマ細胞 (NS-1) を融合して得られたハイブリドーマ
Isotype	IgG ₁

Application notes ELISA
Recommended use

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Storage 2~8°C 3ヶ月間
-20°C 3年間
(凍結融解の繰り返しは避けてください)

References



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619