# RayBio<sup>®</sup> Human GATA-1 Transcription Factor Activity Assay Kit

Catalog #: TFEH-GATA1

User Manual Mar 13, 2017



3607 Parkway Lane, Suite 200 Norcross, GA 30092

Tel: 1-888-494-8555 (Toll Free) or 770-729-2992, Fax: 770-206-2393 Web: www.RayBiotech.com, Email: info@raybiotech.com



# RayBio® Human GATA-1 TF-Activity Assay Kit Protocol TABLE OF CONTENTS

| l.    | Introduction                  | 2   |
|-------|-------------------------------|-----|
| II.   | How It Works                  | 4   |
| III.  | Reagents                      | 5   |
| IV.   | Storage                       | 5   |
| V.    | Additional Materials Required | 6   |
| VI.   | Reagent Preparation           | 6   |
| VII.  | Assay Procedure               | 8   |
| VIII. | Assay Procedure Summary       | 10  |
| IX.   | Typical Data                  | 11  |
| Χ.    | Troubleshooting Guide         | .13 |

Please read the entire manual carefully before starting your experiment.

#### I. INTRODUCTION

GATA transcription factors have a characteristic feature that they all bind to the DNA consensus sequence (A/T)GATA(A/G) of target gene promoters by two adjacent zinc fingers (Cys-X2-Cys-X17-Cys-X2-Cys). These transcription factors regulate differentiation, growth and survival of a wide range of cell types. The GATA family consists of six transcription factors, GATA1 to GATA6, between which the DNA-binding regions, two zinc finger domains are highly conserved. The C-terminal "C-finger" and its adjacent basic region are necessary and sufficient for GATA to bind its cognate sequence, The N-terminal "N-finger" can also bind DNA independently but has a preference for GATC core motifs. Both fingers participate in binding to the palindromic GATA motif ATCWGATA (W = A/T), resulting in markedly increased affinity.

Based on the expression profiles and function in different cell types, GATA family members are divided into two subfamilies. GATA1, GATA2, and GATA3 belong to the hematopoietic subfamily, since they are expressed mainly in the hematopoietic system. The activation of GATA1 is essential for erythroid and megakaryocytic development. GATA2 plays an essential role in regulating the transcription of genes involved in the development and proliferation of hematopoietic. GATA3 is an important regulator of T cell development, including Th2 and regulatory T cells, and plays an important role in endothelial cell biology. The nonhematopoietic subfamily is composed of GATA4, GATA5, and GATA6, which are predominantly expressed in endodermand mesoderm-derived tissues, including intestine, lung, and heart. GATA-4 and -6 regulate expression of several cardiac-specific genes, and during embryonic development, GATA-4 is essential for proper cardiac morphogenesis. GATA-5 performs a unique temporally and spatially restricted function in the embryonic heart and lung. Moreover, GATA-5 may play an important role in the smooth muscle cell diversity.


Accurate monitoring of the level of activated GATA-1 in cells, tissues or animal models is required for both basic science research

investigating signal transduction pathways and applications such as drug development, and simple, speedy and high-throughput methods are needed for this purpose.

Traditionally, western blot to detect the expression of GATA-1, electrophoretic mobility shift assay (EMSA) to detect the DNA binding capacity of GATA-1, or transfection of reporter genes such as luciferase and  $\beta$ -galactosidase with GATA-1 binding sites in culture cells are used in evaluation of GATA-1 reactivity. However, these methods are time consuming, laborious, and sometimes require the use of radioactivity.

The RayBio®® GATA-1 TF-Activity Assay (Transcription Factor-Activity Assay) kit is a non-radioactive transcription factor activity assay with an ELISA format. It offers an easy, speedy, sensitive and high-throughput method to detect the activation of transcription factors. In 96-well plates, double stranded oligonucleotides containing GATA-1 binding sequence have been coated. These oligonucleotides specifically capture the active GATA-1 contained in whole cell lysate or nuclear extracts after a short incubation. Subsequently, the primary antibody against GATA-1 recognizes the GATA-1-DNA complex in each well, and a HRP-conjugated secondary antibody is then used for detection. After washing away any unbound antibody, signal can be obtained easily through a colorimetric assay with a spectrophotometric plate reader at 450 nm. The specificity of the reaction between active GATA-1 and the DNA probe is additionally stringent because of the establishment of specific competitive DNA and non-specific competitive DNA probes in this reaction system.

# **II. HOW IT WORKS**



### **III. STORAGE**

Upon receipt, the positive control should be removed and stored at -20° or -80°C. The remainder of the kit can be stored for up to 6 months at 2-8°C from the date of shipment. Opened Microplate Wells or reagents may be stored for up to 1 month at 2-8°C. Return unused wells to the pouch containing desiccant pack, reseal along entire edge.

Note: The kit can be used within one year if the whole kit is stored at -20°C upon receipt. Avoid repeated freeze-thaw cycles.

### **IV. REAGENTS**

| Component                            | Description   | Size            |
|--------------------------------------|---|-----------------|
| GATA-1 DNA Probe<br>Microplate       | 96 wells (12 strips X 8 wells) coated with GATA-1 probes                              | 1 plate         |
| DNA Binding Buffer                   | 5X concentrated Buffer  | 4 ml            |
| Positive Control                     | Cell nuclear extracts   | 1 vial (20 µl)  |
| Specific Competitor DNA<br>Probe     | Free DNA probes that compete with the coated probes by binding with activated GATA-1. | 1 vial          |
| Non-specific Competitor<br>DNA Probe | Free DNA probes with mutations of the coated DNA probe. Cannot bind activated GATA-1. | 1 vial          |
| Assay Reagent                        | 1X solution   | 1 vial (200 µl) |
| DTT                                  | 300 mM DTT  | 1 vial (200 µl) |
| Wash Buffer Concentrate (20X)        | 20X concentrated solution   | 25 ml           |
| GATA-1 Primary Antibody              | Anti- GATA-1 antibody   | 1 vial          |
| HRP-conjugated Secondary Antibody    | Anti-IgG HRP conjugated antibody  | 1 vial          |
| Antibody Diluent Buffer              | Buffer solution for diluting primary and secondary antibodies                         | 25 ml           |
| TMB One-Step Substrate Reagent       | 3,3,5,5'-tetramethylbenzidine (TMB) in buffer solution                                | 12 ml           |
| Stop Solution                        | 0.2 M sulfuric acid   | 8 ml            |

### V. ADDITIONAL MATERIALS REQUIRED

- 1 Microplate reader capable of measuring absorbance at 450 nm.
- 2 Precision pipettes to deliver 1 µl to 1 ml volumes.
- 3 Adjustable 1-25 ml pipettes for reagent preparation.
- 4 100 ml and 1 liter graduated cylinders.
- 5 Absorbent paper.
- 6 Distilled or deionized water.
- 7 Tubes to prepare positive or sample mixtures.

#### VI. REAGENT PREPARATION

- 1. Preparation of samples:
  - Prepare nuclear extraction or whole lysate containing targeted protein GATA-1 from cell culture or tissue. We recommend using the RayBiotech Nuclear Extraction Kit (Cat#: NE-50) to isolate nuclear proteins for subsequent use in this transcription factor assay.
- 2. Preparation of transcription factor binding reaction system: Bring all reagents to room temperature (18 - 25°C) before use. Thaw the positive control and samples and keep them on ice before adding into wells. Prepare 100 μl transcription factor binding reaction system for each well with 5 x TF-Activity Assay DNA Binding Buffer, TF-Activity Assay Reagent, DTT, Specific Competitor DNA Probe, Non-specific Competitor DNA Probe, and Positive Control or samples containing targeted proteins. Typical examples are shown in the table below.

# <u>Note:</u>

Each reaction may be prepared in a labeled microfuge tube or directly in the coated plate well. If the reaction system is prepared directly in the coated plate wells, please add the reagents sequentially as shown in the table to get the best results.

|   | REACTION  |   |   |   |   |
|---|---|---|---|---|---|
| COMPONENT                                     | Positive control  | Sample  | Specific competitor   | Non-Specific competitor   | Blank   |
| 5x TF-Activity<br>Assay DNA<br>Binding Buffer | 20 μΙ   | 20 μΙ   | 20 μl   | 20 μΙ   | 20 μΙ   |
| TF-Activity Assay Reagent                     | 1.5 µl  | 1.5 μl  | 1.5 µl  | 1.5 µl  | 1.5 μΙ  |
| DTT   | 1 μΙ  |
| Specific Competitor                           | 1   | -   | 10 μl   | -   |   |
| Non-specific<br>Competitor                    | 1   | -   | -   | 10 µl   | -   |
| Control/Sample containing proteins            | 5 μΙ  | * µl  | * µl  | * µl  | -   |
| Total volume                                  | bring final<br>volume to<br><b>100µl</b> with<br>deionized<br>water | bring final<br>volume to<br><b>100μl</b> with<br>deionized<br>water | bring final<br>volume to<br><b>100μl</b> with<br>deionized<br>water | bring final<br>volume to<br><b>100µl</b> with<br>deionized<br>water | bring final<br>volume to<br><b>100µl</b> with<br>deionized<br>water |

<sup>\*</sup> Please note that the amount of total protein containing the target protein to be used in this test can be optimized and must be determined by the investigator.

# Briefly spin down the TF-Activity Assay GATA-1 Primary Antibody vial. Add 60 µl of Antibody Diluent Buffer into the vial to prepare a primary antibody concentrate. Pipette up and down to mix gently (the concentrate can be stored at 4°C for 5 days). The primary antibody concentrate should then be diluted 100-fold

with the Antibody Concentrate should then be diluted 100-fold with the Antibody Diluent Buffer and used in step 4 of Part VII Assay Procedure.

4. Preparation of secondary antibody: Briefly spin down the TF-Activity Assay HRP-conjugated Secondary Antibody vial before use. Add 100 μl of Antibody Diluent Buffer into the vial to prepare a detection antibody

3. Preparation of primary antibody:

concentrate. Pipette up and down to mix gently (the concentrate can be stored at 4°C for 5 days). The detection antibody concentrate should then be diluted 100-fold with the Antibody Diluent Buffer and used in step 6 of Part VII Assay Procedure.

# 5. Preparation of 1x Wash Buffer:

Dilute 25 ml of the 20x Wash Buffer Concentrate into deionized or distilled water to yield 500 ml of 1x Wash Buffer. If the Wash Buffer Concentrate (20x) contains visible crystals, warm to room temperature and mix gently until dissolved.

### Note:

All reagents containing protein (positive control, samples) should be kept on ice to maintain protein stability.

If the reaction system is prepared directly in the coated plate wells, please add the reagents sequentially as shown in the table to get the best results.

To observe the specificity of the DNA binding activity, the amount of protein used in wells of sample, specific competitor and non-specific competitor must be the same.

A positive control should be included every time to confirm correct operation of experiment, however it is not necessary to run specific competitor and non-specific competitor for each sample and every time.

### VII. ASSAY PROCEDURE:

1. Bring the 96-well plate to room temperature (18 - 25°C) before use. If the whole plate will not be used in this assay, place remaining wells back to 2 to 8°C or -20°C. It is recommended that all positive control and samples be run at least in duplicate.

- 2. Add 100 µl of each prepared transcription factor binding reaction system (see Reagent Preparation step 2) including positive control, specific competitor, non-specific competitor and sample into appropriate wells. Cover wells and incubate for 2 hours at room temperature or overnight at 4°C with gentle shaking.
- 3. Discard the solution and wash 4 times by filling each well with 300  $\mu$ l of 1x Wash Buffer (Reagent Preparation step 5) using a multi-channel pipette or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
- 4. Add 100 μl of prepared TF-Activity Assay GATA-1 Primary Antibody (Reagent Preparation step 3) to each well. Incubate for 1 hour at room temperature with gentle shaking.
- 5. Discard the solution. Repeat the wash as in step 3.
- 6. Add 100 μl of prepared TF-Activity Assay HRP-conjugated Secondary Antibody (see Reagent Preparation step 4) to each well. Incubate for 1 hour at room temperature with gentle shaking.
- 7. Discard the solution. Wash as directed in step 3.
- 8. Add 100  $\mu$ l of TMB One-Step Substrate Reagent (Item H) to each well. Incubate for 30 minutes at room temperature in the dark with gentle shaking.
- 9. Add 50  $\mu$ l of Stop Solution (Item I) to each well. Read at 450 nm immediately.

# **VII. ASSAY PROCEDURE SUMMARY**

1. Prepare all reagents, samples, and controls as instructed.

 $\prod$ 

2. Add 100 μl reaction mixture to each well. Incubate 2 hours at room temperature or overnight at 4°C.

 $\prod$ 

3. Add 100 µl prepared primary antibody to each well. Incubate 1 hour at room temperature.

 $\int$ 

4. Add 100 μl prepared secondary antibody. Incubate 1 hour at room temperature.

 $\int$ 

5. Add 100  $\mu$ l TMB One-Step Substrate Reagent to each well. Incubate 30 minutes at room temperature.

 $\prod$ 

6. Add 50 μl Stop Solution to each well. Read at 450 nm immediately.

## VIII. TYPICAL DATA




Figure 1: Transcription factor activity assay of GATA-1 from nuclear extracts of K562 cells or HeLa cells. In K562 cells, activated GATA-1 is translocated into the nucleus where it binds with its corresponding DNA. A. Western-blot result of GATA-1 from cytoplasmic and nuclear fractions. B. Transcription factor activity assay of GATA-1 from

nuclear fractions with the RayBio® GATA-1 TF-Activity Assay Kit (cat # TFEH- GATA-1).




Fig. 2: Transcription factor activity assay of GATA-1 from nuclear extracts of K562 cells or HeLa cells with the specific competitor or non-specific competitor. The result shows specific binding of GATA-1 to the GATA conserved binding site detected by using the RayBio<sup>®</sup> GATA-1 TF-Activity Assay Kit (cat # TFEH- GATA1).

# TROUBLESHOOTING GUIDE

| Problem            | Cause   | Solution  |
|--------------------|---|---|
| 1. Low signal      | 1.Too brief incubation times  | Ensure sufficient incubation time; assay procedure step 2 change to overnight   |
|                    | <ol> <li>Missed key reagent,<br/>inadequate reagent<br/>volumes or improper<br/>dilution</li> </ol> | <ol> <li>Check all reagents         have been added         and check pipettes         to ensure correct         preparation</li> </ol> |
|                    | Not enough targeted protein per well  | <ol> <li>Check positive<br/>control wells and<br/>increase the amount<br/>of sample.</li> </ol>   |
|                    | <ol> <li>Inadequate<br/>development in<br/>colorimetric assay</li> </ol>                            | <ol> <li>Ensure correct<br/>developing buffer<br/>and enough time<br/>used</li> </ol>   |
| 2. Large CV        | 1. Inaccurate pipetting   | 1. Check pipettes   |
|                    | Wells cross     contamination   | Be careful when<br>preparing samples<br>between wells   |
| 3. High background | Plate is insufficiently     washed  | Review the manual for proper wash. If using a plate washer, check that all ports are unobstructed                                       |
|                    | Contaminated wash     Buffer  | Make fresh wash     buffer  |
|                    | Incorrect antibody     dilution   | Check antibody dilutions  |

# RayBio® TF-Activity Assay kits:

Choose TF-Activity Assay kits with more targets for human, mouse, rat and a variety of other species. Visit <a href="https://www.raybiotech.com">www.raybiotech.com</a> for the complete list.

This product is for research use only

