

Mouse IL-6 ELISA Kit

Catalog Number: CEA-M254

Assay Tests: 96 tests

For Research Use Only. Not For Use in Diagnostic or Therapeutic Procedures

CEA-M254-EN03

IMPORTANT: Please carefully read this user guide before performing your experiment.

Product information

This kit is specifically designed for the accurate quantitation of mouse IL-6 from cell culture supernates, serum

and plasma.

The principle of this assay employs a quantitative sandwich enzyme immunoassay approach. Initially, a microplate

is coated with a capture antibody. Then, samples and biotinylated capture antibody are added to the wells. After

the removal of any unbound materials through washing, streptavidin-HRP (SA-HRP) conjugate is added to the

wells. Streptavidin has a very high affinity for biotin, so it binds to the biotinylated capture antibody that is already

bound to the target antigen. After washing, a substrate specific to HRP is added to the wells. HRP catalyzes a

reaction that converts the substrate into a detectable signal, often a color change or luminescence, depending

on the substrate used. This enzymatic reaction amplifies the signal, allowing for higher sensitivity in detecting the

target analyte. The intensity of the signal is measured using a spectrophotometer.

.

NOTE:

1. This kit is for research use only and is not for use in diagnostic or therapeutic applications.

2. Please do not use the kit after the expiration date indicated on the kit label.

3. Do not mix or substitute reagents with those from other lots or sources.

Manufactured and distributed by

ACRODiagnostics Inc.

TEL: 010-67855298-8147

Distributed by:

ACROBiosystems Inc.

US & Canada TEL: +1 800-810-0816

Asia & Pacific TEL: +86 400-682-252

Contents

The kit contains sufficient reagents for 96 wells.

Catalog	Contents	Amount
CEA254-C01	Pre-coated Anti-mouse-IL-6 Antibody Microplate	1 plate
CEA254-C02	Mouse IL-6 Standard	6 ng×2
CEA254-C03	Biotin-Anti-mouse-IL-6 Antibody	2.0 μg
CEA254-C04	Biotin-Antibody Dilution Buffer	8 mL
CEA254-C05	Streptavidin-HRP Con. Solution	500 μL
CEA254-C06	Streptavidin-HRP Dilution Buffer	15 mL
CEA254-C07	20× Washing Buffer	50 mL
CEA254-C08	Sample Dilution Buffer	15 mL×2
CEA254-C09	Substrate Solution	12 mL
CEA254-C10	Stop Solution	6 mL

NOTE: Bubbles in microplate wells do not affect the experiment and require no action. Proceed with the experimental procedures and methods described below.

Storage

Keep the unopened kit stored at 2-8 °C. Avoid using the kit beyond its expiration date. For opened kit and reconstituted reagents, with the exception of the two contents listed in following table, others can be stored for up to 30 days at 2-8 °C.

Contents	Storage conditions
Pre-coated Anti-mouse-IL-6 Antibody Microplate	Return unused wells to the foil pouch, reseal along entire edge. May be stored for up to 1 month at 2-8°C.
Mouse IL-6 Standard	Aliquot and store for up to 1 month at -70°C in a freezer. Avoid repeated freeze-thaw cycles.
Biotin-Anti-mouse-IL-6 Antibody	Aliquot and store for up to 1 month at -70°C in a freezer. Avoid repeated freeze-thaw cycles.

NOTE: Streptavidin-HRP Con. Solution and Substrate Solution should avoid light.


IMPORTANT: Bring all reagents to room temperature before use. If crystals have formed in buffer solution, place the buffer solution in an 37°C incubator until the crystals have completely dissolved and bring the solution back to room temperature before use.

Required materials not supplied.

Instrument	Microplate reader capable of measuring absorbance at 450 nm
Reagents	Deionized, ultrapure or distilled water
	50 mL and 500 mL graduated cylinders
Consumables	Pipettes and pipette tips
	Tubes to prepare standard dilutions.

Workflow

Analyte: Mouse IL-6



NOTE: Incubation temperature is 18 $^{\circ}$ C-25 $^{\circ}$ C

Prepare the working buffers and standard dilutions.

Prepare the reconstituted Biotin-Anti-mouse-IL-6 Antibody

Add 100 $\,\mu$ L ultrapure water to the provided lyophilized product (CEA254-C03) , dissolve at room temperature for 15-30 minutes, and mix by gently pipetting. The concentration of reconstituted Biotin-Anti-mouse-IL-6 Antibody Solution is 20.0 $\,\mu$ g/mL.

NOTE: Avoiding vigorous shaking or vortexing. The reconstituted solution should be stored at -70°C. The freeze-thaw cycle should not exceed 1 time.

Prepare the working buffers.

- 1. 1×Washing Buffer: Dilute 50 mL 20×Washing Buffer with deionized or distilled water to 1000 mL.
- 2. Biotin-Anti-mouse-IL-6 Antibody Solution: Add 60 μ L of the reconstituted Biotin-Anti-mouse-IL-6 Antibody stock Solution to 6 mL Biotin-Antibody Dilution Buffer, thoroughly mix. The solution was freshly prepared just before use.
- 3. Mouse IL-6 Streptavidin-HRP Solution: Add 400 μ L of Streptavidin-HRP Con. Solution to 12 mL of Streptavidin-HRP Dilution Buffer, thoroughly mix. The solution was freshly prepared just before use.

Prepare the reconstituted standard.

Add 1 mL ultrapure water to the provided lyophilized product (CEA254-C02) , dissolve at room temperature for 15-30 minutes, and mix by gently pipetting. The concentration of reconstituted mouse IL-6 Standard is 6 ng /mL.

NOTE: Avoiding vigorous shaking or vortexing. The reconstituted solution should be stored at -70°C. The freeze-thaw cycle should not exceed 1 time.

Prepare the standard serial dilutions.

- 1. Label 7 tubes, one for each standard point: Std.-1, Std.-2, Std.-3, Std.-4, Std.-5, Std.-6, Std.-7.
- 2. Add 100 μ L of the liquid from reconstituted mouse IL-6 Standard and 900 μ L of Sample Dilution Buffer to tube Std.-1, thoroughly mix (Std.-1 = 600 pg/mL).
- 3. Prepare 1:1 serial dilutions for the standard curve as follows: Add 500 μL of Sample Dilution Buffer to each tube (Std.-2, Std.-3, Std.-4, Std.-5, Std.-6, Std.-7).
- 4. Transfer 500 μ L of liquid from Std.-1 to the tube Std.-2, and thoroughly mix (Std.-2 = 300 pg/mL).
- 5. Continue to transfer 500 μ L of liquid from previous dilution tube to the next dilution tube until add liquid to tube Std.-7.
- 6. Sample Dilution Buffer serves as zero standard (blank).

PROCEDURE OF ASSAY

- 1. Add 50 μ L of mouse IL-6 Standard, sample, or control to wells. Seal the plate with microplate sealing film. Incubate at room temperature (18-25 °C) for **1 hour**.
- 2. Aspirate each well and add 300 μ L of 1×Washing Buffer to each well, gently tap the plate for **1 minute**. Remove any remaining Washing Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels. Repeat the wash process four times for a total of five washes.
- 3. Add 50 μL Biotin-Anti-mouse-IL-6 Antibody Solution to each well, Seal the plate with microplate sealing film. Incubate at room temperature (18-25 °C) for **1 hour.**
- 4. Repeat step 2.
- 5. Add 100 μ L of Mouse IL-6 Streptavidin-HRP Solution to each well. Seal the plate with microplate sealing film. Incubate at room temperature (18-25 °C) for 30 minutes, avoid light.
- 6. Repeat step 2.
- 7. Add 100 μ L of Substrate Solution to each well. Seal the plate with microplate sealing film and incubate at room temperature (18-25 °C) for **15 minutes, avoid light**.
- Add 50 μL of Stop Solution to each well. Tap the plate gently to ensure thorough mixing.
 Note: the color in the wells should change from blue to yellow.
- 9. Read the absorbance at 450nm and 630nm using Microplate reader within 10minutes.

 *Note: To reduce the background noise, subtract the readings at 630nm from the readings at 450nm.

CALCULATION OF RESULTS

- 1. Compute the average of the duplicated readings for every standard, control, and sample. Then, subtract the average optical density (O.D.) of the zero standard(blank).
- 2. Establish a standard curve by processing the data using computer software capable of executing a **four-parameter logistic (4-PL)** curve fitting.
- 3. Normal range of Standard curve: $R^2 \ge 0.9900$.
- 4. If the OD value of the sample to be tested is higher than the highest standard, the sample shall be diluted with dilution buffer and assay repeated.

Typical data

Note: For each experiment, a standard curve needs to be set for each microplate, and the specific OD value may vary depending on different laboratories, testers, or equipment. The following example data is for reference only. The sample concentration was calculated based on the results of the standard curve.

louse IL-6 Standard (pg/mL)	OD _{450nm-630nm}	R ² =0.9998
600	2.555	37
300	1.588	ns ity
150	0.896	a de la companya de l
75	0.461	Optical Density
37.5	0.272	° lee
18.75	0.179	0 200 400 600 800 Conc.[pg/mL]
Blank	0.089	


PERFORMANCE CHARACTERISTICS

1. Sensitivity

The minimum detectable concentration (MDC) of Mouse IL-6 is typically less than 10 pg/mL. The MDC was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

2. Linearity

Three samples (Serum) spiked with high concentrations of mouse IL-6 were serially diluted with dilution buffer to produce samples with values within the dynamic range of the assay and then assayed. The average recovery of mouse IL-6 for serum samples is 97.3%.



3. Intra-Assay Precision

Ten replicates of each of 3 samples containing different mouse IL-6 concentrations were tested in one assay. Acceptable criteria: CV < 10%.

Sample Concentration (pg/mL)	Mean (pg /mL)	SD	Numbers	CV (%)
500	515.03	34.42	10	6.7
125	129.15	6.65	10	5.2
25	26.49	1.57	10	5.9

4. Inter-Assay Precision

Five samples containing different concentrations of mouse IL-6 were tested in independent assays. Acceptable criteria: CV<15%.

Sample Concentration (pg/mL)	Mean (pg/mL)	SD	Numbers	CV (%)
500	502.09	36.85	9	7.3
125	127.92	9.24	9	7.2
25	26.69	2.67	9	10.0

5. Recovery


Recombinant mouse IL-6 was spiked into 3 mouse serum samples, and then analyzed. The average recovery of mouse IL-6 for serum samples is 88.7%.

Sample ID	Conc Measured (pg/mL)	Conc Added (pg/mL)	Conc Recovered (pg/mL)	Recovery (%)
	285.44	300	283.55	94.5
1	130.05	150	128.17	85.5
	68.07	75	66.18	88.3
	2.10	-		
	266.18	300	264.29	88.1

	125.28	150	123.39	82.3
2	66.48	75	64.59	86.1
	4.18	-		
	302.46	300	300.57	100.2
3	137.88	150	135.99	90.7
3	73.70	75	71.82	95.8
	6.36	-		

6. Hook Effect

Not be affected by the concentration of mouse IL-6 up to 5 ng/mL.



7. Sample Values

The mice in the experimental group were injected with 15 μ g or 30 μ g of LPS and bled 4 hours after injection. In contrast, those in the control group received an equivalent volume of normal saline (NS). Serum samples were collected from all mice, and IL-6 concentrations were quantitatively measured using standardized immunoassays.

Condition	Mean (ng/mL)
Injected with NS	0.011
Injected with 15 μg LPS	484.89
Injected with 30 µg LPS	320.61

8. Specificity

No cross-reactivity was observed when this kit was used to analyze the following recombinant cytokines at up to 1 μ g/mL.

Mouse	IL-1β、IL-2、TNF-α、IFN-γ
Human	IL-6

9. CALIBRATION

This immunoassay is calibrated against highly purified recombinant mouse IL-6 produced at ACROBiosystems. The NIBSC/WHO International Standard for mouse IL-6 (93/730), which was intended as a potency standard, was evaluated in this kit. NIBSC/WHO (93/730) approximate value (U/mL) = $0.25 \times ACRO$ Mouse IL-6(pg/mL)

TROUBLESHOOTING GUIDE

Problem	Cause	Solution	
Poor standard curve	* Inaccurate pipetting	* Check pipettes	
Large CV	* Inaccurate pipetting* Air bubbles in wells	* Check pipettes* Remove bubbles in wells	
High background	* Plate is insufficiently washed* Contaminated wash buffer	* Review the manual for proper wash. * Make fresh wash buffer	
Very low readings across the plate * Incorrect wavelengths * Insufficient development time		* Check filters/reader * Increase development time	
Samples are reading too high, but standard curve looks fine * Samples contain cytokine levels above assay range		* Dilute samples and run again	
Drift	* Interrupted assay set-up * Reagents not at room temperature	* Assay set-up should be continuous - have all standards and samples prepared appropriately before commencement of the assay * Ensure that all reagents are at room temperature before pipetting into the wells unless otherwise instructed in the antibody inserts	