



Human Tumor Necrosis Factor Alpha (TNF-α) ELISA Kit (Enzyme-Linked Immunosorbent Assay)

Catalog Number: CRS-A002

Pack Size: 96 tests

IMPORTANT: Please carefully read this manual before performing your experiment.

For Research Use Only. Not For Use in Diagnostic or Therapeutic Procedures



ClinMax A002-EN.05

INTENDED USE

The kit is developed for quantitative detection of TNF- α in human serum and cell culture supernates. It is intended for research use only (RUO).

BACKGROUND

Tumor Necrosis Factor alpha (TNF- α), also known as cachectin, is the prototypic ligand of the TNF superfamily. It is produced primarily by activated macrophages, also secreted by other cells such as CD4+lymphocyte, NK cell, Neutrophils, mast cells, eosinophils and neurons. It is a pleiotropic molecule and play a central role on inflammation, immune system development, apoptosis and lipid metabolism. Its overexpression related to a series of pathological state including Cachexia, septic shock, and autoimmune disorders.

PRINCIPLE OF THE ASSAY

This assay kit is used to measure the levels of human Tumor Necrosis Factor alpha (TNF- α) by employing a standard sandwich-ELISA format. The micro-plate in the kit has been pre-coated with Anti-TNF- α Antibody. Firstly, add the standard samples provided in kit and your samples to the plate, incubate and wash the wells. Then add the Biotin-Anti-TNF- α Antibody to the plate and form Antibody-antigen-biotinylated antibody complex, incubate and wash the wells. Next add Streptavidin-HRP to the plate, incubate and wash the wells. At last, load the substrate into the wells and monitor solution color from blue to yellow. The reaction is stopped by the addition of a stop solution and the intensity of the absorbance can be measured at 450nm and 630nm. The OD Value reflects the amount of TNF- α bound.

PRECAUTIONS

- 1. This kit is for research use only and is not for use in diagnostic or therapeutic applications.
- 2. The kit is suitable for cell supernatant, serum and plasma samples.
- 3. Do not use reagents past their expiration date.
- 4. Do not mix or substitute reagents with those from other kits or other lot number kits.
- 5. If samples generate values higher than the highest standard, dilute the samples with the appropriate

_

US and Canada: Tel: +1 800-810-0816 Web: http://www.acrobiosystems.com





4002-EN 05

calibrator diluent and repeat the assay. If cell supernatant samples need step dilution, except for the final dilution with diluent, other intermediate dilutions can be in cell culture medium.

- 6. Differences in test results can be caused by a variety of factors, including laboratory operator, pipette usage, plate washing technique, reaction time or temperature, and kit storage.
- 7. This kit is designed to remove or reduce some endogenous interference factors in biological samples, and not all possible influencing factors have been removed.

MATERIALS PROVIDED

Table 1. Materials provided

Catalog	Size Forma		Format	Stor	orage	
Cilling	Components	(96 tests)	1 01 mat	Unopened	Opened	
CRS002-C01	Pre-coated Anti-TNF-α Antibody Microplate	1 plate	Solid	2-8°C	2-8°C	
CRS002-C02	Human TNF-α Standard	20 μg	Powder	2-8°C	-70°C	
CRS002-C03	Biotin-Anti- TNF-α Antibody	20 μg	Powder	2-8°C	-70°C	
CRS002-C04	Streptavidin-HRP	50 μL	Liquid	2-8°C, avoid light	2-8°C, avoid light	
CRS002-C05	10xWashing Buffer	50 mL	Liquid	2-8°C	2-8°C	
CRS002-C06	2xDilution Buffer	50 mL	Liquid	2-8°C	2-8°C	
CRS002-C07	Substrate Solution	12 mL	Liquid	2-8°C, avoid light	2-8°C, avoid light	
CRS002-C08	Stop Solution	7 mL	Liquid	2-8°C	2-8°C	

2

US and Canada: Tel: +1 800-810-0816 Web: http://www.acrobiosystems.com





KIT STORAGE AND EXPIRATION DATE

- 1. The unopened kit is stable for 18 months from the date of manufacture if stored at 2°C to 8°C.
- 2. The opened kit should be stored per Table 1. The shelf life is 30 days from the date of opening.

Note: a. Do not use reagents past their expiration date.

b. Find the expiration date on the outside packaging.

REAGENTS/EQUIPMENT NEEDED BUT NOT SUPPLIED

Single or multi-channel micropipettes and pipette tips: need to meet $10 \mu L$, $300 \mu L$, $1000 \mu L$ injection; requirements;

37° C Incubator;

Single or dual wavelength microplate reader with 450nm and 630nm filter;

Tubes: 1.5mL,10mL;

Timer;

Reagent bottle;

Deionized or distilled water.

REAGENT PREPARATION

Bring all reagents and samples to room temperature (20°C-25°C) before use. If crystals have formed in buffer solution, place the sample in an 37°C incubator until the crystals have completely dissolved and bring the solution back to room temperature before use.

According to Table 2, prepare the provided lyophilized product into a storage solution with ultrapure water, dissolve at room temperature for 15 to 30 minutes, and mix by gently pipetting, avoiding vigorous shaking or vortexing. The reconstituted storage solution should be stored at -70° C. It is recommended that the number of freezing and thawing should not exceed 1 time, and the size of the aliquot should not be less than $10 \mu g$.

Table 2. Preparation method

Catalog	Components	Size (96 T)	Storage solution concentration.	Reconstituted water Vol.	
CRS002-C02	Human TNF-α Standard	20 μg	100 μg/mL	200 μL	
CRS002-C03	Biotin-Anti-TNF-α Antibody	20 μg	100 μg/mL	200 μL	

3

Web: http://www.acrobiosystems.com

US and Canada: Tel: +1 800-810-0816





RECOMMENDED SAMPLE PREPARATION

1. Working Solution Preparation

1.1 Preparation of 1×Washing Buffer:

Dilute 50 mL 10×Washing Buffer with ultrapure water/deionized water to 500 mL.

1.2 Preparation of 1×Dilution Buffer:

Dilute 50 mL 2×Dilution Buffer with 1×Washing Buffer to 100 mL.

1.3 Preparation of Biotin-Anti-TNF-α Antibody working fluid:

Dilute Biotin-Anti-TNF-α Antibody with 1×Dilution Buffer to 0.5 µg/mL. Please prepare it for one-time use only.

1.4 Preparation of Streptavidin-HRP working fluid:

Dilute Streptavidin-HRP at 1:2000 with 1×Dilution Buffer. The prepared working fluid should avoid light. Please prepare it for one-time use only.

- 1.5 Sample preparation
- a. If the sample to be tested is the serum or plasma, dilute test sample at 1:2 with 1×Dilution Buffer. The volume ratio of sample to diluent is 1:1.
- b. If the sample to be tested is the cell supernatant, dilute test sample at 1:2 with 1×Dilution Buffer. The volume ratio of sample to diluent is 1:1.

2. Preparation of Standard curve

The concentration of the reconstituted human TNF-α Standard (CRS002-C02) is 100 μg/mL, prepare (Std.-0) by diluting 5 μL the reconstituted human TNF-α Standard into 495 μL Sample Dilution Buffer, mix gently well. Then prepare Std.-1' by diluting 5 μL Std.-0 into 495 μL Sample Dilution Buffer. At last, prepare the highest concentration of standard curve, Std.-1 (1250 pg/mL), by diluting 75 μL Std.-1' into 525 μL Sample Dilution Buffer. Prepare 1:1 serial dilutions for the standard curve as follows: Pipette 300 μL of Sample Dilution Buffer into each tube. Make sure to mix well every time. Sample Dilution Buffer serves as blank.


4

US and Canada: Tel: +1 800-810-0816 Web: http://www.acrobiosystems.com





A002-EN.05



3. Add Samples

Add 100 μ L Calibrator and samples to each well. For blank Control wells, please add 100 μ L Dilution Buffer.

Note: It is recommended to set doeble holes for samples and standard curves to be tested.

4. Incubation

Seal the plate with microplate sealing film and incubate at room temperature for 1 hour.

5. Washing

Remove the remaining solution by aspiration, add 300 µL of 1×Washing Buffer to each well, soak for 10s, remove any remaining 1×Washing Buffer: by aspirating or decanting, invert the plate and blot it against paper towels. Repeat the wash step above for three times.

6. Add Biotin-Anti-TNF-α Antibody

For all wells, add 100 μ L Biotin- Anti-TNF- α Antibody (0.5 μ g/mL) working solution. Please prepare it for one-time use only.

7. Incubation

Seal the plate with microplate sealing film and incubate at room temperature for 1 hour.

Tel: +86 400-682-2521

8. Washing

Repeat step 5.

9. Add Streptavidin-HRP

Asia and Pacific:

5

US and Canada: Tel: +1 800-810-0816

E-mail: order@acrobiosystems.com

Web: http://www.acrobiosystems.com





For all wells, add $100 \mu L$ Streptavidin-HRP (dilute at 1:2000) working solution. Please prepare it for one-time use only, avoid light.

10. Incubation

Seal the plate with microplate sealing film and incubate at room temperature for 30 min.

11. Washing

Repeat step 5.

12. Substrate Reaction

Add 100 µL Substrate Solution to each well. Seal the plate with microplate sealing film and incubate at room temperature for 20 min, avoid light.

13. Termination

Add 50 µL Stop Solution to each well, and tap the plate gently to allow thorough mixing.

Note: The color in the wells should change from blue to yellow.

14. Data Recording

Read the absorbance at 450 nm and 630 nm using UV/Vis microplate spectrophotometer within 10 minutes

Note: To reduce the background noise, subtract the readings at 630nm from the readings at 450nm.

CALCULATION OF RESULTS

- 1. Calculate the mean absorbance for each standard, control and sample and subtract average zero standard optical density (O.D.).
- 2. The standard curve is plotted with the standard concentration as x-axis and the calibrated absorbance value as y-axis. Four parameters logistic are used to draw the standard curve and calculate the sample concentration.
- 3. Normal range of Standard curve: R²≥0.9900.
- 4. Detection range: 19.5 pg/mL-1250 pg/mL. If the OD value of the sample to be tested is higher than 1250 pg/mL, the sample shall be diluted with dilution buffer and assay repeated. If the OD value of the sample to be tested is lower than 19.5 pg/mL, the sample should be reported.

6

E-mail: order@acrobiosystems.com

US and Canada: Tel: +1 800-810-0816 Web: http://www.acrobiosystems.com

Asia and Pacific: Tel: +86 400-682-2521





QUICK GUILD



7

US and Canada: Tel: +1 800-810-0816

Asia and Pacific: Tel: +86 400-682-2521

E-mail: order@acrobiosystems.com


Web: http://www.acrobiosystems.com





TYPICAL DATA

The following data is for reference only. The sample concentration was calculated based on the results of the standard curve.



SENSITIVITY

The minimum detectable concentration of human TNF-α is 2.321 pg/mL. The minimum detectable concentration was determined by adding twice standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

PRECISION

1. Intra-assay Precision

Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision.

2. Precision

Three samples of known concentration were tested on 3 different plates, 10 replicates in each plate to assess inter-assay precision.

Intra-assay Precision	Inter-assay Precision
-----------------------	-----------------------

8

US and Canada:

Tel: +1 800-810-0816

Web: http://www.acrobiosystems.com
E-mail: order@acrobiosystems.com

Asia and Pacific: Tel:

Tel: +86 400-682-2521





A	n	N	2	_1	F)	N	•	n	5	
$^{-}$	v	v	'	-,	Ŀ.	LN	٠,	U	J	

						1100 <u>=</u> E100
Sample	1	2	3	1	2	3
n	20	20	20	30	30	30
Mean (pg/mL)	346.018	172.273	89.184	359.774	184.266	97.260
SD	8.042	11.376	2.782	17.914	10.386	7.461
CV (%)	2.3	6.6	3.1	5.0	5.6	7.7

Note: The example data is for reference only.

RECOVERY

Three parts of blank serum were added with different concentrations of human TNF- α , and the serum without human TNF- α was used as background to calculate the recovery rate. The Range of the recovery rate is 90%-111%, and the average recovery is 101%.

Sample Type	Average% Recovery	Range		
Serum(n=5)	101	90-111%		

LINEARITY

To assess the linearity of the assay, samples spiked with high concentrations of human TNF- α were serially diluted with calibrator diluent to produce samples with values within the dynamic range of the assay.

		Cell culture medium	Serum	Citrate plasma
1.2	Average Recovery (%)	109	105	102
1:2	Range (%)	102-116	96-114	96-111
1.4	Average Recovery (%)	110	101	96
1:4	Range (%)	105-114	95-113	91-102
1.0	Average Recovery (%)	97	89	105
1:8	Range (%)	95-101	86-96	96-112

Note: The example data is for reference only.

SPECIFICITY

This assay recognizes natural and recombinant human TNF- α . No cross-reactivity was observed when this kit was used to analyze the following recombinant cytokines.

9

US and Canada: Tel: +1 800-810-0816

Asia and Pacific: Tel: +86 400-682-2521

E-mail: order@acrobiosystems.com

Web: http://www.acrobiosystems.com





Human					
IL-1β	IL-10				
IL-2	IL-15				
IL-4	IL-21				
IL-7	GM-CSF				
IL-6	IFN-gamma				

SAMPLE VALUES

a. 80 human healthy sample were evaluated for the presence of human TNF-α in this assay. No medical histories were available for the donors used in this study.

Sample	n	Concentration	Detection percentage (%)	Detection average Concentration
Sample	n	(pg/mL)	Detection percentage (%)	(pg/mL)
Serum	80	n.d 295.711	35	62.534

Note: n.d. means that the concentration value cannot be detected which is lower than the sensitivity. The example data is for reference only.

b. In this assay, the sample is the cell supernatant of PBMC proliferation stimulated by CD28/OKT3 antibody.

Sampla	n	Concentration	Detection percentage	Detection average Concentration
Sample	n	(pg/mL)	(%)	(pg/mL)
cell supernatant	2	102.613-230.675	100	166.644

Note: The example data is for reference only.

CALIBRATION

This immunoassay is calibrated against a highly purified *E. coli*-expressed recombinant human TNF- α (17/232). Reference Reagent is calibrated by NIBSC/WHO in June 2020.

NIBSC/WHO (17/232) approximate value (U/mL) = $0.110 \times \text{Human TNF-}\alpha \text{ value (pg/mL)}$

PLATE LAYOUT

10

US and Canada: Tel: +1 800-810-0816 Web: http://www.acrobiosystems.com
Asia and Pacific: Tel: +86 400-682-2521 E-mail: order@acrobiosystems.com





	1	2	3	4	5	6	7	8	9	10	11	12
А	Std1	Std1	(···)							()		
В	Std2	Std2	$\left(\cdots \right)$	()	()	()		()	()	()		
С	Std3	Std3	$\left(\begin{array}{c} \cdots \end{array} \right)$	()		()	()	()	()	()		
D	Std4	Std4		()	()	()		()	()	()		
E	Std5	Std5	$\left(\cdots \right)$	())		()	()	())	
F	Std6	Std6	()	()	()	()		()		()		
G	Std7	Std7	()	()	()	()	()	()	()	()	()	
н	Blank	Blank	()	()	()	···)			()	()		

Note: Blank is a Blank Dilution Buffer hole.

11

US and Canada: Tel: +1 800-810-0816 Web: http://www.acrobiosystems.com





TROUBLESHOOTING GUIDE

Problem	Cause	Solution
Poor standard curve	* Inaccurate pipetting	* Check pipettes
Large CV	* Inaccurate pipetting	* Check pipettes
	* Air bubbles in wells	* Remove bubbles in wells
High background	* Plate is insufficiently	* Review the manual for proper wash.
	washed	* Make fresh wash buffer
	* Contaminated wash buffer	
Very low readings	* Incorrect wavelengths	*Check
across the plate	* Insufficient development	filters/reader
	time	* Increase
		development
		time
Samples are reading too	* Samples contain cytokine	* Dilute samples and run again
high, but standard	levels above assay range	
curve looks fine		
Drift	* Interrupted assay set-up	* Assay set-up should be continuous - have all
	* Reagents not at room	standards and samples prepared appropriately
	temperature	before commencement of theassay
		* Ensure that all reagents are at room temperature
		before pipetting into the wells unless otherwise
		instructed in the antibody inserts

12

Web: http://www.acrobiosystems.com US and Canada: Tel: +1 800-810-0816 E-mail: order@acrobiosystems.com

Asia and Pacific: Tel: +86 400-682-2521