## Biotinylated SARS-CoV-2 Spike RBD Protein, His, Avitag™ (BF.7&BA.4.6/Omicron) (MALS verified)

Catalog # SPD-C82E1





#### **Synonym**

Spike,S protein RBD,Spike glycoprotein Receptor-binding domain,S glycoprotein RBD,Spike protein RBD

#### Source

Biotinylated SARS-CoV-2 Spike RBD, His,Avitag (BF.7&BA.4.6/Omicron) (SPD-C82E1) is expressed from human 293 cells (HEK293). It contains AA Arg 319 - Lys 537 (Accession # QHD43416.1 (G339D, R346T, S371F, S373P, S375F, T376A, D405N, R408S, K417N, N440K, L452R, S477N, T478K, E484A, F486V, Q498R, N501Y, Y505H)). The spike mutations are identified on the SARS-CoV-2 Omicron variant (Pango lineage: BF.7&BA.4.6). Predicted N-terminus: Arg 319

## **Molecular Characterization**

This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (Avitag<sup>TM</sup>).

The protein has a calculated MW of 28.3 kDa. The protein migrates as 34-36 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

#### Labeling

Biotinylation of this product is performed using Avitag<sup>TM</sup> technology. Briefly, the single lysine residue in the Avitag is enzymatically labeled with biotin.

#### **Protein Ratio**

Passed as determined by the HABA assay / binding ELISA.

## **Purity**

>95% as determined by SDS-PAGE.

>90% as determined by SEC-MALS.

### **Formulation**

Lyophilized from  $0.22~\mu m$  filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

#### Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

#### Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:


- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

#### **SDS-PAGE**



Biotinylated SARS-CoV-2 Spike RBD, His, Avitag (BF.7&BA.4.6/Omicron) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%.

#### SEC-MALS



The purity of Biotinylated SARS-CoV-2 Spike RBD, His, Avitag (BF.7&BA.4.6/Omicron) (Cat. No. SPD-C82E1) is more than 90% and the molecular weight of this protein is around 34-44 kDa verified by SEC-MALS. Report

# **Bioactivity-ELISA**




## Biotinylated SARS-CoV-2 Spike RBD Protein, His,Avitag™ (BF.7&BA.4.6/Omicron) (MALS verified)

Catalog # SPD-C82E1






# Biotinylated SARS-CoV-2 Spike RBD, His, Avitag (BF.7&BA.4.6/Omicron) ELISA 0.5 μg of Human ACE2, Fc Tag per well



Biotinylated SARS-CoV-2 Spike RBD, His, Avitag (BF.7&BA.4.6/Omicron) Conc. (ng/mL)

Immobilized Human ACE2, Fc Tag (Cat. No. AC2-H5257) at 5  $\mu$ g/mL (100  $\mu$ L/well) can bind Biotinylated SARS-CoV-2 Spike RBD, His,Avitag<sup>TM</sup> (BF.7&BA.4.6/Omicron) (Cat. No. SPD-C82E1) with a linear range of 0.2-6 ng/mL (QC tested).


Biotinylated SARS-CoV-2 Spike RBD, His, Avitag (BF.7&BA.4.6/Omicron) ELISA 0.1  $\mu$ g of Biotinylated SARS-CoV-2 Spike RBD, His, Avitag (BF.7&BA.4.6/Omicron) per well



Anti-SARS-CoV-2 Spike RBD Antibody, Chimeric mAb, Human IgG1 (AM130) Conc. (ng/mL)

Immobilized Biotinylated SARS-CoV-2 Spike RBD, His,Avitag<sup>TM</sup> (BF.7&BA.4.6/Omicron) (Cat. No. SPD-C82E1) at 1 μg/mL (100 μL/well) on streptavidin (Cat. No. STN-N5116) precoated (0.5 μg/well) plate can bind Anti-SARS-CoV-2 Spike RBD Antibody, Chimeric mAb, Human IgG1 (Cat. No. S1N-M130) with a linear range of 0.1-1 ng/mL (Routinely tested).

# Biotinylated SARS-CoV-2 Spike RBD, His, Avitag (BF.7&BA.4.6/Omicron) ELISA 0.1 μg of Biotinylated SARS-CoV-2 Spike RBD, His, Avitag (BF.7&BA.4.6/Omicron) per well



Anti-SARS-CoV-2 Spike RBD Broadly Neutralizing Antibody, Human IgG1 Conc. (ng/mL)

Immobilized Biotinylated SARS-CoV-2 Spike RBD, His,Avitag<sup>TM</sup> (BF.7&BA.4.6/Omicron) (Cat. No. SPD-C82E1) at 1 μg/mL (100 μL/well) on streptavidin (Cat. No. STN-N5116) precoated (0.5 μg/well) plate can bind Anti-SARS-CoV-2 Spike RBD Broadly Neutralizing Antibody, Human IgG1 (Cat. No. SPD-M265) with a linear range of 0.1-1 ng/mL (Routinely tested).

## Background

It's been reported that Coronavirus can infect the human respiratory epithelial cells through interaction with the human ACE2 receptor. The spike protein is a large type I transmembrane protein containing two subunits, S1 and S2. S1 mainly contains a receptor binding domain (RBD), which is responsible for recognizing the cell surface receptor. S2 contains basic elements needed for the membrane fusion. The S protein plays key parts in the induction of neutralizing-antibody and T-cell responses, as well as protective immunity.

