

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK  
Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: [info@abbexa.com](mailto:info@abbexa.com)

## **Snailase**

Catalogue No.: abx082089

Snailase is an enzyme.

**Target:** Snailase

**Conjugation:** Unconjugated

**Directions for use:** After processing the sulfhydryl compounds, suspend the yeast cells in isotonic sorbitol solution (pH 5.8-7.2). Add 30-40 mg enzyme per 1 g cells, incubate at 37 °C for 1 hour to dissolve the cytoderm. The efficiency of fragmentation is greater than 90%.

**Note:** This product is for research use only.