

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK  
Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

## **Coomassie Brilliant Blue Fast Staining Solution**

Catalogue No.: abx090652

Based on Coomassie G250 dye, Coomassie Brilliant Blue Fast Staining Solution can be used for fast and sensitive staining proteins in electrophoresis gel in SDS-PAGE or native PAGE, or used for detecting the remaining proteins on PAGE after Western transfer membrane. Instead of common toxic methanol and irritant acetic acid, this solution is prepared by a unique formula, it is non-toxic and non-irritant. This solution can detect band at 100ng in electrophoresis gel within 1 hour or shorter time, however, the routine method need more than 3 hours.

**Tested Applications:**

WB

**Note:**

This product is for research use only.