

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK
Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

PCR Stimulant

Catalogue No.: abx098001

PCR Stimulant can be used to increase sensitivity and specificity for GC/AT-rich template or complex template. It is especially suitable for Pfu enzymes. The stock concentration is 5x, and the working concentration can be varied between 0.5x and 2.5x.

Target: PCR Stimulant

Storage: Store at -20°C. Stable for 12 months from date of receipt.

Note: This product is for research use only.