

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK
Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

PCR SuperMix for PAGE (+dye)

Catalogue No.: abx098004

PCR SuperMix for PAGE is a ready-to-use mixture of DNA Polymerase for PAGE, dNTPs and optimized buffer. The supermix for PAGE is provided at 2x concentration and used at 1x concentration by adding template, primers and H₂O.

Target: PCR SuperMix for PAGE (+dye)

Storage: Store at -20°C. Stable for 12 months from date of receipt.

Note: This product is for research use only.