

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK
Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

cDNA Synthesis SuperMix for qPCR

Catalogue No.: abx098019

First-Strand cDNA Synthesis SuperMix for qPCR provides all the necessary components for cDNA synthesis from total RNA or mRNA. It is provided at 5x concentration and used at 1x concentration by adding gDNA remover, RNA and H₂O. The resulting cDNA are suitable for qPCR and not for regular PCR.

Target: First-Strand cDNA Synthesis SuperMix for qPCR (One-Step gDNA Removal)

Storage: Store at -20°C. Stable for 12 months from date of receipt.

Note: This product is for research use only. This product is shipped with dry ice.