

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK
Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

miRNA First-Strand cDNA Synthesis SuperMix

Catalogue No.: abx098020

miRNA First-Strand cDNA Synthesis SuperMix provides all the necessary components for cDNA synthesis from miRNA. High efficient poly(A) tail addition and first-strand cDNA synthesis are performed miRNA RT Enzyme Mix (containing tailing enzyme and RT enzyme) and 2xTS miRNA Reaction Mix. This kit permits simple procedure to realize reverse transcription of miRNA-containing samples like total RNA or small RNAs. During the procedure, other RNAs like mRNA and snRNA can also be reverse transcribed.

Target: miRNA First-Strand cDNA Synthesis SuperMix

Storage: Store at -20°C. Stable for 12 months from date of receipt.

Note: This product is for research use only. This product is shipped with dry ice.