

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK
Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

One-Step RT-PCR SuperMix (+dye) (4 kb)

Catalogue No.: abx098027

One-Step RT-PCR combines the first-strand cDNA synthesis with PCR in the same tube to simplify reaction setup and reduce the possibility of contamination. Only gene-specific primers can be used for One-Step RT-PCR. Amplified fragment up to 4 kb.

Target: One-Step RT-PCR SuperMix (+dye) (4 kb)

Storage: Store at -20°C. Stable for 12 months from date of receipt.

Note: This product is for research use only. This product is shipped with dry ice.