

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK
Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

Ni-NTA Resin

Catalogue No.: abx098108

Ni-NTA Resin is used for one-step purification of Histagged proteins. The His-tagged proteins bind to Ni²⁺ cations, which are immobilized on the Ni-NTA resin by 4 metal-chelating sites. After unbound proteins are washed away, the target proteins are recovered by gradient elution. It is suitable for both native and denatured protein purification.

Target: Ni-NTA Resin

Storage: Store at -20°C. Stable for 12 months from date of receipt.

Note: This product is for research use only.