

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK
Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

RNase-Free water

Catalogue No.: abx098145

RNase-free Water is prepared with deionized water, which is treated with 0.01% DEPC, followed by high-temperature autoclave. It is suitable for a wide range of RNA-related molecular biology applications.

Target: RNase-Free water

Storage: Store at -20°C. Stable for 12 months from date of receipt.

Note: This product is for research use only.