

## **DATASHEET**

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Tax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

Catalogue No.:abx098432

Haptoglobin Assay Kit is an Assay Kit by Immunoturbidimetry Method.

Based on the bottle type, product model can be classified into Hitachi 7170, Hitachi 7060, Hitachi 7020, Beckman, Toshiba, Mindray, DuPont, Innova, Siemens, KHB, Abbott, Universal types, etc.

This product is used to determine haptoglobin content in human serum.

The HP in human serum meets the specific antibody (goat anti - human HP) in the liquid buffer, forming antigen - antibody complex, and the turbidity formed is proportional to the HP concentration in the linear range. HP concentration in the sample can be calculated by comparing the calibrator result with same processing.

R1: Detection Reagent 1	R2: Detection Reagent 2
Phosphate Buffer → 20 mmol/l	Goat-anti human HP antibody → Appropriate
PEG, Preservatives → Apprpriate	