

## **DATASHEET**

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK
Telephone: +44 (0) 1229 75595a - Fax: +44 (0) 1327 755951 - E-Mail: info@abbexa.com

Catalogue No.:abx098460

Transferrin Assay Kit is an Assay Kit by Immunoturbidimetry Method.

Based on the bottle type, product model can be classified in to Hitachi 7170, Hitachi 7060, Hitachi 7020, Beckman, Toshiba, Mindray, DuPont, Innova, Siemens, KHB, Abbott, Universal types, etc

This product is used to determine transferrin content in human serum.

The TF in human serum m eets the specific antibody (goat anti - human TF) in the liquid buffer, forming antigen - antibody complex, and the turbidity formed is proportional to the TF concentration in the linear range. TF concentration in the sample can be calculated by comparing the calibrator result with same processing.

R1: Detection Reagent 1	R2: Detection Reagent 2
Phosphate Buffer → 10 mmol/L	Goat Anti-Human TF Antibody → Appropriate
PEG, Preservatives → Appropriate	