

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK
Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

PCR SuperMix

Catalogue No.: abx098887

PCR SuperMix is a ready-to-use mixture of DNA polymerase, dNTPs, and optimized buffer. The SuperMix is provided at 2x concentration and can be used at 1x concentration by adding template, primers and H₂O. PCR SuperMix offers 108-fold fidelity as compared to DNA Polymerase. The Extension rate is about 2-6 kb/min. PCR products can be directly cloned into Blunt vectors. Amplification of genomic DNA fragment up to 15 kb. Amplification of plasmid DNA fragment up to 20 kb.

Target: PCR SuperMix

Storage: Store at -20 °C for two years.

Note: This product is for research use only.