

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK
Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

RNA Beads

Catalogue No.: abx298002

RNA beads can easily extract RNA from mastermix of reaction such as rRNA depletion, DNase I digestion, in vitro transcription, labeled RNA and synthetic RNA. The resulting RNA product is suitable for RNA library construction, RT-PCR, qRT-PCR, chip analysis, Northern blot RNAi, or other downstream applications.

Target:	RNA Beads
Storage:	Store at 2-8 °C for one year.
Note:	This product is for research use only.