

## **DATASHEET**

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

## Interferon Gamma (IFNg) Antibody Pair

Catalogue No.:abx370093

Interferon Gamma (IFNg) Antibody Pair for use in Sandwich ELISA assay development. This antibody pair contains:

- Interferon Gamma (IFNg) capture antibody,
- Interferon Gamma (IFNg) biotin-conjugated detection antibody,
- Interferon Gamma (IFNg) standard.

It is recommended to use this antibody pair with abx098958 Antibody Pair Support Kit (Sandwich Method).

Target:	Interferon Gamma (IFNg)
Reactivity:	Human
Tested Applications:	ELISA
Recommended dilutions:	Dilute the Capture Antibody with Coating Buffer.  Dilute the biotin-conjugated Detection Antibody with Detection Antibody Diluent.  Optimal dilutions/concentrations should be determined by the end user.
Form:	Liquid (Capture Antibody and Detection Antibody)
Reconstitution:	Reconstitute the standard with Standard Diluent. The volume, and therefore standard concentration, should be determined by the end user.
Storage:	Store at 2 to 8 °C for up to one month. Aliquot and store at -80 °C for up to one year. Avoid repeated freeze/thaw cycles. All solutions should be made fresh before the experiment.
Standard Form:	Lyophilized
ELISA Type:	Sandwich
Capture Antibody Conjugation:	Unconjugated
Detection Antibody Conjugation:	Biotin
Buffer:	The capture and detection antibody both contain 0.1% sodium azide.



**DATASHEET** 

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

## Directions for use:

Bring all components to room temperature (18-25°C) and briefly spin or centrifuge the vials before use. Working solutions should be prepared and used immediately.

## Recommended Procedure:

- 1. Dilute the Capture Antibody to working concentration using Coating Buffer. Immediately coat the 96-well plate with diluted Capture Antibody (100  $\mu$ l per well). Seal the plate and incubate at 4 °C overnight or at 37 °C for 2 hours
- 2. Aspirate the wells and wash with Wash Buffer (350 µl per well) and allow to soak for 1-2 min. Remove the liquid by inverting and tapping the plate on to absorbent paper.
- 3. Block the plate with Blocking Buffer (200 µl per well) at 37 °C for 1.5 hours.
- 4. Repeat the aspiration/wash process in Step 2.
- 5. Add 100  $\mu$ l of standards or sample into the appropriate wells. Cover with a plate sealer and incubate at 37 °C for 1 hour.
- 6. Repeat the aspiration/wash process in Step 2.
- 7. Add appropriately diluted biotin-conjugated Detection Antibody (100 µl per well). Cover the plate with a new plate sealer and incubate at 37 °C for 1 hour.
- 8. Repeat the aspiration/wash process in Step 2.
- 9. Add appropriately diluted Streptavidin HRP (100  $\mu$ l per well). Cover the plate with a new plate sealer and incubate at 37 °C for 30 min.
- 10. Repeat the aspiration/wash process in Step 2.
- 11. Add Substrate Solution (90  $\mu$ l per well). Cover the plate with a new plate sealer and incubate at 37 °C for 10-20 min. Keep the plate in the dark and avoid exposure to light.
- 12. Add Stop Solution (50 µl per well). Tap the side of the plate to ensure thorough mixing.
- 13. Measure the absorbance immediately using a microplate reader set at 450 nm.

This product is for research use only.

Note: