

DATASHEET

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

Mouse Phosphoribosylformylglycinamidine Synthase (PFAS) ELISA Kit

Catalogue No.:abx390234



Mouse Phosphoribosylformylglycinamidine Synthase (PFAS) ELISA Kit is an ELISA kit against Phosphoribosylformylglycinamidine Synthase (PFAS).

Target: Phosphoribosylformylglycinamidine Synthase (PFAS)

Reactivity: Mouse

Tested Applications: ELISA

Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Test Range: 78 pg/ml - 5000 pg/ml

Sensitivity: < 33 pg/ml

Validity: The validity for this kit is 6 months.

Storage: Store at 4 °C upon receipt.

Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% within

the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested

that the whole assay is performed by the same user throughout.

Swiss Prot: Q5SUR0

GeneID: 237823

Gene Symbol: PFAS

Ensembl: ENSMUSG00000020899

Standard Form: Lyophilized



DATASHEET

Abbexa Ltd, Innovation Centre, Cambridge Science Park, Cambridge, CB4 0EY, UK Telephone: +44 (0) 1223 755950 - Fax: +44 (0) 1223 755951 - E-Mail: info@abbexa.com

ELISA Detection: Colorimetric

ELISA Type: Sandwich

ELISA Data: Quantitative

Note: This product is for research use only.

The range and sensitivity is subject to change. Please contact us for the latest product information. For accurate results, sample concentrations must be diluted to mid-range of the kit. If you require a

Please note that our ELISA and CLIA kits are optimised for detection of native samples, rather than

specific range, please contact us in advance or write your request in your order comments.

recombinant proteins. We are unable to guarantee detection of recombinant proteins, as they may

have different sequences or tertiary structures to the native protein.