



Technical support: support@abbkine.com

Website: https://www.abbkine.com

## Anti-GAPDH Mouse Monoclonal Antibody (2B5), AbFluor™ 555 Conjugated

Cat #: A01020A555 Size: 100µl

## **Product Information**

	Product Name: Anti-GAPDH Mouse Monoclonal Antibody (2B5), AbFluor™ 555 Conjugated		
	Applications: IHC, IF		Isotype: Mouse IgG1
	Reactivity: Human, Mouse, Rat, Monkey, Dog, Chicken, Rabbit, Pig, Sheep, Insect, Yeast, Hamster		
REF	Catalog Number: A01020A555	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
Î	Storage: Store at -20°C protected from light.	Λ	Note: Contain sodium azide.
1	Avoid repeated freeze / thaw cycles.	<b>\(\cdot \cdot \cd</b>	

**Background:** Glyceraldehyde 3-phosphate dehydrogenase (abbreviated as GAPDH or less commonly as G3PDH) is an enzyme of ~37kDa that catalyzes the sixth step of glycolysis and thus serves to break down glucose for energy and carbon molecules. In addition to this long established metabolic function, GAPDH has recently been implicated in several non-metabolic processes, including transcription activation, initiation of apoptosis ER to Golgi vesicle shuttling, and fast axonal, or axoplasmic transport.

<u>Application Notes</u>: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IHC (1:200), IF (1:200).

Storage Buffer: Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

**Storage Instructions:** Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing. Store in the dark.






Fig. AbFluor™ 555 (λEX/λEm: 555/565 nm) is a bright and photostable fluorescent dye which is super alternative to Cy3, Tetramethylrhodamine, Dylight 549 and Alexa Fluor 555.

<u>Note:</u> The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

