

EME1 Antibody (Center) Blocking Peptide
Synthetic peptide
Catalog # BP6288c**Specification****EME1 Antibody (Center) Blocking Peptide -
Product Information**Primary Accession [Q96AY2](#)**EME1 Antibody (Center) Blocking Peptide -
Additional Information**

Gene ID 146956

Other NamesCrossover junction endonuclease EME1,
3122-, MMS4 homolog, hMMS4, EME1,
MMS4**Target/Specificity**

The synthetic peptide sequence used to generate the antibody [AP6288c](/product/products/AP6288c) was selected from the Center region of human EME1. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.

Format

Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.

Storage

Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.

Precautions

This product is for research use only. Not for use in diagnostic or therapeutic procedures.

**EME1 Antibody (Center) Blocking Peptide -
Protein Information**

Name EME1

**EME1 Antibody (Center) Blocking Peptide -
Background**

EME1 and MUS81 (MIM 606591) form an endonuclease complex that cleaves branched DNA structures, especially those arising during stalled DNA replication (Abraham et al., 2003 [PubMed 14609959]).[supplied by OMIM]

**EME1 Antibody (Center) Blocking Peptide -
References**

Taylor,E.R., Proc. Natl. Acad. Sci. U.S.A. 105 (10), 3757-3762 (2008)Olsen,J.V., Cell 127 (3), 635-648 (2006)

Synonyms MMS4**Function**

Interacts with MUS81 to form a DNA structure-specific endonuclease with substrate preference for branched DNA structures with a 5'-end at the branch nick. Typical substrates include 3'-flap structures, replication forks and nicked Holliday junctions. May be required in mitosis for the processing of stalled or collapsed replication forks.

Cellular Location

Nucleus, nucleolus. Note=Recruited to regions of DNA damage in S-phase cells

EME1 Antibody (Center) Blocking Peptide - Protocols

Provided below are standard protocols that you may find useful for product applications.

- [Blocking Peptides](#)