

#### **ECE1 Antibody (Center) Blocking Peptide**

Synthetic peptide Catalog # BP6855d

#### **Specification**

ECE1 Antibody (Center) Blocking Peptide - Product Information

Primary Accession P42892

ECE1 Antibody (Center) Blocking Peptide - Additional Information

**Gene ID 1889** 

#### **Other Names**

Endothelin-converting enzyme 1, ECE-1, ECE1

#### **Target/Specificity**

The synthetic peptide sequence used to generate the antibody <a href=/products/AP6855d>AP6855d</a> was selected from the Center region of human ECE1. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.

#### Format

Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.

#### Storage

Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.

#### **Precautions**

This product is for research use only. Not for use in diagnostic or therapeutic procedures.

ECE1 Antibody (Center) Blocking Peptide - Protein Information

Name ECE1

**Function** 

## ECE1 Antibody (Center) Blocking Peptide - Background

ECE-1 is involved in proteolytic processing of endothelin precursors to biologically active peptides.

## ECE1 Antibody (Center) Blocking Peptide - References

Cottrell,G.S., et.al., J. Biol. Chem. 284 (33), 22411-22425 (2009)





Converts big endothelin-1 to endothelin-1.

# **Cellular Location**Cell membrane; Single-pass type II membrane protein

#### **Tissue Location**

All isoforms are expressed in umbilical vein endothelial cells, polynuclear neutrophils, fibroblasts, atrium cardiomyocytes and ventricles. Isoforms A, B and C are also expressed in placenta, lung, heart, adrenal gland and phaeochromocytoma; isoforms A and C in liver, testis and small intestine; isoform B, C and D in endothelial cells and umbilical vein smooth muscle cells; isoforms C and D in saphenous vein cells, and isoform C in kidney

### ECE1 Antibody (Center) Blocking Peptide - Protocols

Provided below are standard protocols that you may find useful for product applications.

• Blocking Peptides