

Anti-Integrin α V β 6 (ITGAV & ITGB6) Reference Antibody (STX-100)
Recombinant Antibody
Catalog # APR10321**Specification**

Anti-Integrin α V β 6 (ITGAV & ITGB6) Reference Antibody (STX-100) - Product Information

Application	FC, Kinetics, Animal Model
Primary Accession	P06756 & P18564
Reactivity	Human
Clonality	Monoclonal
Isotype	IgG1
Calculated MW	145.12 KDa

Anti-Integrin α V β 6 (ITGAV & ITGB6) Reference Antibody (STX-100) - Additional Information**Target/Specificity**Integrin α V β 6 (ITGAV & ITGB6)**Endotoxin**< 0.001EU/ μ g,determined by LAL method.**Conjugation**

Unconjugated

Expression system

CHO Cell

Format


Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.

Anti-Integrin α V β 6 (ITGAV & ITGB6) Reference Antibody (STX-100) - Protein Information**Anti-Integrin α V β 6 (ITGAV & ITGB6) Reference Antibody (STX-100) - Protocols**


Provided below are standard protocols that you may find useful for product applications.

- [Western Blot](#)
- [Blocking Peptides](#)
- [Dot Blot](#)
- [Immunohistochemistry](#)
- [Immunofluorescence](#)
- [Immunoprecipitation](#)
- [Flow Cytometry](#)
- [Cell Culture](#)

Anti-Integrin α V β 6 (ITGAV & ITGB6) Reference Antibody (STX-100) - Images



Anti-Integrin α V β 6 (ITGAV & ITGB6) Reference Antibody (STX-100) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-Integrin α V β 6 (ITGAV & ITGB6) Reference Antibody (STX-100) is more than 95%, determined by SEC-HPLC.