

**Anti-TNFRSF12A / TWEAKR / CD266 Reference Antibody (enavatuzumab)
Recombinant Antibody
Catalog # APR10399****Specification**

Anti-TNFRSF12A / TWEAKR / CD266 Reference Antibody (enavatuzumab) - Product Information

Application	FC, Kinetics, Animal Model
Primary Accession	Q9NP84
Reactivity	Human
Clonality	Monoclonal
Isotype	IgG1
Calculated MW	150 KDa

Anti-TNFRSF12A / TWEAKR / CD266 Reference Antibody (enavatuzumab) - Additional Information**Target/Specificity**

TNFRSF12A / TWEAKR / CD266

Endotoxin

< 0.001EU/ µg,determined by LAL method.

Conjugation

Unconjugated

Expression system

CHO Cell

Format

Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.

Anti-TNFRSF12A / TWEAKR / CD266 Reference Antibody (enavatuzumab) - Protein Information**Name** TNFRSF12A**Synonyms** FN14**Function**

Receptor for TNFSF12/TWEAK. Weak inducer of apoptosis in some cell types. Promotes angiogenesis and the proliferation of endothelial cells. May modulate cellular adhesion to matrix proteins.

Cellular Location

Membrane; Single-pass type I membrane protein.

Tissue Location


Highly expressed in heart, placenta and kidney. Intermediate expression in lung, skeletal muscle and pancreas

Anti-TNFRSF12A / TWEAKR / CD266 Reference Antibody (enavatuzumab) - Protocols


Provided below are standard protocols that you may find useful for product applications.

- [Western Blot](#)
- [Blocking Peptides](#)
- [Dot Blot](#)
- [Immunohistochemistry](#)
- [Immunofluorescence](#)
- [Immunoprecipitation](#)
- [Flow Cytometry](#)
- [Cell Culture](#)

Anti-TNFRSF12A / TWEAKR / CD266 Reference Antibody (enavatuzumab) - Images



Anti-TNFRSF12A / TWEAKR / CD266 Reference Antibody (enavatuzumab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%



The purity of Anti-TNFRSF12A / TWEAKR / CD266 Reference Antibody (enavatuzumab) is more than 98.7%, determined by SEC-HPLC.