

**Anti-ANGPTL8 Reference Antibody (Regeneron patent anti-ANGPTL8)
Recombinant Antibody
Catalog # APR10784****Specification**

Anti-ANGPTL8 Reference Antibody (Regeneron patent anti-ANGPTL8) - Product Information

Application	FC, Kinetics, Animal Model
Primary Accession	Q6UXH0
Reactivity	Human
Clonality	Monoclonal
Isotype	IgG2SA
Calculated MW	150 KDa

Anti-ANGPTL8 Reference Antibody (Regeneron patent anti-ANGPTL8) - Additional Information**Target/Specificity**
ANGPTL8**Endotoxin**
< 0.001EU/ µg,determined by LAL method.**Conjugation**
Unconjugated**Expression system**
CHO Cell**Format**
Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.**Anti-ANGPTL8 Reference Antibody (Regeneron patent anti-ANGPTL8) - Protein Information****Name** ANGPTL8 ([HGNC:24933](#))**Function**
Hormone that acts as a blood lipid regulator by regulating serum triglyceride levels (PubMed:22569073, PubMed:22809513, PubMed:23150577). May be involved in the metabolic transition between fasting and refeeding: required to direct fatty acids to adipose tissue for storage in the fed state (By similarity).**Cellular Location**
Secreted.

Tissue Location


Predominantly expressed in liver. Also expressed in adipose tissues.

Anti-ANGPTL8 Reference Antibody (Regeneron patent anti-ANGPTL8) - Protocols


Provided below are standard protocols that you may find useful for product applications.

- [Western Blot](#)
- [Blocking Peptides](#)
- [Dot Blot](#)
- [Immunohistochemistry](#)
- [Immunofluorescence](#)
- [Immunoprecipitation](#)
- [Flow Cytometry](#)
- [Cell Culture](#)

Anti-ANGPTL8 Reference Antibody (Regeneron patent anti-ANGPTL8) - Images



Anti-ANGPTL8 Reference Antibody (Regeneron patent anti-ANGPTL8) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-ANGPTL8 Reference Antibody (Regeneron patent anti-ANGPTL8) is more than 95% ,determined by SEC-HPLC.