

**Anti-Complement C5 Reference Antibody (crovalimab)  
Recombinant Antibody  
Catalog # APR10840****Specification**

---

**Anti-Complement C5 Reference Antibody (crovalimab) - Product Information**


|                   |                            |
|-------------------|----------------------------|
| Application       | FC, Kinetics, Animal Model |
| Primary Accession | <a href="#">P01031</a>     |
| Reactivity        | Human                      |
| Clonality         | Monoclonal                 |
| Isotype           | IgG1                       |
| Calculated MW     | 145 KDa                    |

**Anti-Complement C5 Reference Antibody (crovalimab) - Additional Information****Target/Specificity**  
Complement C5**Endotoxin**  
< 0.001EU/ µg,determined by LAL method.**Conjugation**  
Unconjugated**Expression system**  
CHO Cell**Format**  
Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.**Anti-Complement C5 Reference Antibody (crovalimab) - Protein Information****Name** C5**Synonyms** CPAMD4**Function**  
Activation of C5 by a C5 convertase initiates the spontaneous assembly of the late complement components, C5-C9, into the membrane attack complex. C5b has a transient binding site for C6. The C5b-C6 complex is the foundation upon which the lytic complex is assembled.**Cellular Location**  
Secreted.**Anti-Complement C5 Reference Antibody (crovalimab) - Protocols**


Provided below are standard protocols that you may find useful for product applications.

- [Western Blot](#)
- [Blocking Peptides](#)
- [Dot Blot](#)
- [Immunohistochemistry](#)
- [Immunofluorescence](#)
- [Immunoprecipitation](#)
- [Flow Cytometry](#)
- [Cell Culture](#)

#### **Anti-Complement C5 Reference Antibody (crovalimab) - Images**



Anti-Complement C5 Reference Antibody (crovalimab) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-Complement C5 Reference Antibody (crovalimab) is more than 95% ,determined by SEC-HPLC.