# ABclomal www.abclonal.com

# PER2 Rabbit pAb

Catalog No.: A24439

# **Basic Information**

#### **Observed MW**

150-160kDa

#### **Calculated MW**

137kDa

#### Category

Polyclonal Antibody

### **Applications**

WB

### **Cross-Reactivity**

Human, Mouse

# **Background**

This gene is a member of the Period family of genes and is expressed in a circadian pattern in the suprachiasmatic nucleus, the primary circadian pacemaker in the mammalian brain. Genes in this family encode components of the circadian rhythms of locomotor activity, metabolism, and behavior. This gene is upregulated by CLOCK/ARNTL heterodimers but then represses this upregulation in a feedback loop using PER/CRY heterodimers to interact with CLOCK/ARNTL. Polymorphisms in this gene may increase the risk of getting certain cancers and have been linked to sleep disorders.

# **Recommended Dilutions**

**WB** 

1:1000 - 1:5000

## **Immunogen Information**

**Gene ID** 8864

**Swiss Prot** 

015055

### **Immunogen**

Synthetic peptide. This information is considered to be commercially sensitive.

### **Synonyms**

FASPS; FASPS1; PER2

### **Contact**

0

www.abclonal.com

### **Product Information**

Source

**Isotype** 

**Purification** 

Rabbit

IgG


Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.

## **Validation Data**




Western blot analysis of various lysates, using PER2 Rabbit pAb (A24439) at 1:2000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (RM00021).

Exposure time: 60s.



Western blot analysis of lysates from Mouse skeletal muscle, using PER2 Rabbit pAb (A24439) at 1:2000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (RM00021).

Exposure time: 60s.