

# Polyclonal Sheep anti Human Oncostatin M/OSM

## ORDERING INFORMATION

|                 |                                                    |
|-----------------|----------------------------------------------------|
| Catalog Number: | hAP-5369                                           |
| Size:           | 50 µg                                              |
| Storage:        | < -20 °C                                           |
| Immunogen:      | E. coli-derived recombinant human Oncostatin M/OSM |
| Ig Type:        | IgG                                                |
| Purification:   | Antigen Affinity Purified                          |
| Applications:   | WB N/B FC                                          |

### **Preparation:**

This antibody was produced from a sheep immunized with purified human Oncostatin M/OSM protein.

### **Formulation and Storage:**

The IgG fraction of immunized serum was purified by antigen affinity chromatography and lyophilized from a 0.2 µm filtered solution in phosphate-buffered saline (PBS). Lyophilized samples are stable for 2 years from date of receipt when stored at -70°C.

### **Reconstitution:**

Reconstitute the antibody with 100 µl sterile PBS and the final concentration is 500µg/ml. Reconstituted antibody can also be aliquoted and stored frozen at < -20 °C for at least for six months without detectable loss of activity. **Avoid repeated freeze-thaw cycles.**

### **Specificity:**

Detects endogenous human Oncostatin M/OSM in direct ELISAs and Western blots. Less than 1% cross reactivity. Cross reactivity to other species has not been tested!

### **Applications:**

|       |     |                 |
|-------|-----|-----------------|
| WB    | Yes | 1 : 1000 - 2000 |
| N / B | Yes |                 |
| FC    | Yes |                 |

*Optimal dilutions should be determined by each laboratory for each application.*

The listed dilutions are for recommendation only and the final conditions should be optimized by the ender users!

**This product is sold for Research Use Only !**