

DATA SHEET

N-acetylneuraminic acid synthase Human Recombinant

Item Number	rAP-1830
Synonyms	Sialic acid synthase, N-acetylneuraminate synthase, N-acetylneuraminate-9-phosphate synthase, N- acetylneuraminic acid phosphate synthase, N-acetylneuraminic acid synthase, NANS, SAS.
Description	NANS Human Recombinant fused with a 20 amino acid His tag at N-terminus produced in E.Coli is a sin- gle, non-glycosylated, polypeptide chain containing 379 amino acids (1-359 a.a.) and having a molecular mass of 42.4kDa. The NANS is purified by proprietary chromatographic techniques.
Uniprot Accesion Number	Q9NR45
Amino Acid Sequence	MGSSHHHHHH SSGLVPRGSH MPLELELCPG RWVGGQHPCF IIAEIGQNHQ GDLDVAKRMI RMAKECGADC AKFQKSELEF KFNRKALERP YTSKHSWGKT YGEHKRHLEF SHDQYRELQR YAEEVGIFFT ASGMDEMAVE FLHELNVPFF KVGSGDTNNF PYLEKTAKKG RPMVISSGMQ SMDTMKQVYQ IVKPLNPNFC FLQCTSAYPL QPEDVNLRVI SEYQKLFPDI PIGYSGHETG IAISVAAVAL GAKVLERHIT LDKTWKGSDH SASLEPGELA ELVRSVRLVE RALGSPTKQL LPCEMACNEK LGKSVVAKVK IPEGTILTMD MLTVKVGEPK GYPPEDIFNL VGKKVLVTVE EDDTIMEELV DNHGKKIKS.
Source	Escherichia Coli.
Physical Appearance and Stability	Sterile Filtered colorless solution. Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at - 20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).Avoid multiple freeze-thaw cycles.
Formulation and Purity	The NANS solution (0.5 mg/ml) contains 20mM Tris-HCl buffer (pH8.0), 0.1M NaCl, 1mM DTT and 10% glycerol. Greater than 95.0% as determined by SDS-PAGE.
Application	
Solubility	
Biological Activity	
Shipping Format and Condition	Lyophilized powder at room temperature.

Optimal dilutions should be determined by each laboratory for each application. The listed dilutions are for recommendation only and the final conditions should be optimized by the ender users! This product is sold for Research Use Only