



## Fibroblast Growth Factor-Basic Human Recombinant

**Item Number** rAP-2203

Prostatropin, HBGH-2, HBGF-2, FGF-b. Synonyms

Description Fibroblast Growth Factor-2 Human Recombinant (FGF-2) produced in E.Coli is a single, non-glycosylated,

polypeptide chain containing 154 amino acids and having a molecular mass of 17.2kDa.The FGF-b is puri-

fied by proprietary chromatographic techniques.

P09038 **Uniprot Accesion Number** 

AAGSITTLPA LPEDGGSGAF PPGHFKDPKR LYCKNGGFFL RIHPDGRVDG VREKSDPHIK **Amino Acid Sequence** 

LQLQAEERGV VSIKGVCANR YLAMKEDGRL LASKCVTDEC FFFERLESNN YNTYRSRKYT SWY-

VALKRTG QYKLGSKTGP GQKAILFLPM SAKS.

Source Escherichia Coli.

**Physical Appearance** 

and Stability

Sterile Filtered White lyophilized (freeze-dried) powder. Lyophilized Fibroblast Growth Factor-2 although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution FGF-b should be stored at 4°C between 2-7 days and for future use below -18°C.For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).Please prevent freeze-thaw cycles.

The protein was lyophilized from a concentrated (1mg/ml) solution in PBS, pH 7.4. Greater than 98.0% as Formulation and Purity

determined by:(a) Analysis by RP-HPLC.(b) Analysis by SDS-PAGE.

**Application** 

Solubility It is recommended to reconstitute the lyophilized Fibroblast Growth Factor Basic in sterile 18MΩ-cm H2O

not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

The ED50, calculated by the dose-dependant proliferation of BAF3 cells expressing FGF receptors **Biological Activity** 

(measured by 3H-thymidine uptake) is <0.5 ng/ml, corresponding to a specific activity of 2,000,000 Units/

**Shipping Format and Condition** Lyophilized powder at room temperature.

Optimal dilutions should be determined by each laboratory for each application. The listed dilutions are for recommendation only and the final conditions should be optimized by the ender users! This product is sold for Research Use Only