

## Pegylated Mouse Leptin Recombinant

<b>Item Number</b>	rAP-2358
<b>Synonyms</b>	OB Protein, Obesity Protein, OBS, Obesity factor.
<b>Description</b>	Leptin Mono-Pegylated Mouse Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 146 amino acids and an additional Ala at N-terminus. Pegylated Mouse Leptin contains PEG 20 kDa at its N-terminus and having a molecular mass of 35.6 kDa as determined by mass spectrom-
<b>Uniprot Accession Number</b>	P41160
<b>Amino Acid Sequence</b>	
<b>Source</b>	Escherichia Coli.
<b>Physical Appearance and Stability</b>	Sterile Filtered White lyophilized (freeze-dried) powder. Lyophilized Leptin although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Leptin should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.
<b>Formulation and Purity</b>	The mouse Leptin was lyophilized from a concentrated (0.65mg/ml) solution containing 0.003mM NaHCO <sub>3</sub> . Greater than 99.0% as determined by: (a) Analysis by Gel-Filtration. (b) Analysis by SDS-PAGE.
<b>Application</b>	
<b>Solubility</b>	It is recommended to reconstitute the lyophilized Leptin in sterile 0.4% NaHCO <sub>3</sub> adjusted to pH-8.5 not less than 100µg/ml, which can then be further diluted to other aqueous solutions.
<b>Biological Activity</b>	Pegylated mouse Leptin is capable of stimulating proliferation of BAF/3 cells stably transfected with the long form of human leptin receptor. Pegylated mouse Leptin in vitro activity is only slightly lower than the non-pegylated antagonist but in vivo it h
<b>Shipping Format and Condition</b>	Lyophilized powder at room temperature.

Optimal dilutions should be determined by each laboratory for each application. The listed dilutions are for recommendation only and the final conditions should be optimized by the ender users! This product is sold for **Research Use Only**