

DnaK (HSP70) E.Coli Recombinant

Item Number	rAP-3373
Synonyms	HSP-70, HSP70, DnaK, , Chaperone protein dnaK, Heat shock protein 70, Heat shock 70 kDa protein, groP, grpF, seg, b0014, JW0013.
Description	Recombinant DnaK produced in E.Coli is a single, non-glycosylated polypeptide chain containing 638 amino acids and having a molecular mass of 69 kDa.
Uniprot Accession Number	P0A6Y8
Amino Acid Sequence	MGKIIIGIDLG TTNSCVAIMD GTTPRVLENA EGDRTTPSII AYTQDGETLV GQPAKRQAVTNPQNTLFAIK RLIGRRFQDE EVQRDVSIMP FKIIAADNGD AWVEVKGQKM APPQISAEVLKMKKTAEDY LGPEVTEAVI TVPAYFNDAQ RQATKDAGRI AGLEVKRIINEPTAAALAYGLDKGTGNRTI AVYDLGGGTDFDISIIEDEV DGEKTFEVL TNGDTHLGGE DFDSRLINYLVEEFKKDQGI DLRNDPLAMQ RLKEAAEKAK IELSSAQQTD VNLPIYTADA TGPKHMNIKVTRAKLESLE DLVNRSEIPL KVALQDAGLS VSDIDDVILV GGQTRMPMVQ KKVAEFFGKEPRKDVNPDEA VAIGAAVQGG VLTGDVKDVL LLDVTPLSLG
Source	Escherichia Coli.
Physical Appearance and Stability	Sterile filtered colorless solution. Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Formulation and Purity	The DnaK protein contains 25mM Tris-HCl buffer (pH 7.5), 100mM NaCl, 5mM DTT and 10% Glycerol. Greater than 95.0% as determined by SDS-PAGE.
Application	
Solubility	
Biological Activity	
Shipping Format and Condition	Lyophilized powder at room temperature.

Optimal dilutions should be determined by each laboratory for each application. The listed dilutions are for recommendation only and the final conditions should be optimized by the ender users! This product is sold for **Research Use Only**