

EnzyChrom™ Glycogen Assay Kit (Cat# E2GN-100)

Quantitative Colorimetric/Fluorimetric Glycogen Determination

DESCRIPTION

GLYCOGEN is a branched polysaccharide of glucose units linked by α -1,4 glycosidic bonds and α -1,6 glycosidic bonds. It is stored primarily in the liver and muscle, and forms an energy reserve that can be quickly mobilized to meet a sudden need for glucose. The most common glycogen metabolism disorder is found in diabetes, in which, due to abnormal amounts of insulin, liver glycogen can be abnormally accumulated or depleted. Genetic glycogen storage diseases have been associated with various inborn errors of metabolism caused by deficiencies of enzymes necessary for glycogen synthesis or breakdown.

Simple, direct and automation-ready procedures for measuring glycogen concentrations find wide applications in research and drug discovery. BioAssay Systems' glycogen assay uses a single Working Reagent that combines the enzymatic break down of glycogen and the detection of glucose in one step. The color intensity of the reaction product at 570nm or fluorescence intensity at $\lambda_{\text{ex/em}} = 530/585$ nm is directly proportional to the glycogen concentration in the sample. This simple convenient assay is carried out at room temperature and takes only 30 min.

KEY FEATURES

Use as little as 10 μL samples. Linear detection range: 2 to 200 $\mu\text{g/mL}$ glycogen for colorimetric assays and 0.2 to 20 $\mu\text{g/mL}$ for fluorimetric assays.

KIT CONTENTS (100 TESTS IN 96-WELL PLATES)

Assay Buffer: 12 mL **Dye Reagent:** 120 μL
Enzyme A: Dried **Standard:** 50 μL 50 mg/mL
Enzyme B: 120 μL

Storage conditions. The kit is shipped on ice. Store all components at -20°C upon receiving. Shelf life: 6 months after receipt.

Precautions: reagents are for research use only. Normal precautions for laboratory reagents should be exercised while using the reagents. Please refer to Material Safety Data Sheet for detailed information.

PROCEDURES

Reagent Preparation:

Reconstitute Enzyme A by adding 120 μL Assay Buffer to the Enzyme A tube. Make sure Enzyme A is fully dissolved by pipetting up and down. Store reconstituted Enzyme A at -20°C and use within 1 month.

Sample Preparation:

Samples can be prepared according to Murat & Serfaty (*Clin Chem* 20:1576-1577, 1974). Briefly, homogenize tissue/cell sample in 25 mM citrate, pH 4.2, 2.5 g/L NaF on ice. Centrifuge 14,000 g for 5 min to remove debris, and use 10 μL clear supernatant for the assay.

Colorimetric Procedure:

1. Equilibrate all components to room temperature. During experiment, keep thawed enzymes in a refrigerator or on ice.
2. **Standards and samples:** Dilute standard by mixing 5 μL Standard with 1.245 mL dH_2O to give 200 $\mu\text{g/mL}$ standard. Dilute standard in dH_2O as follows.

No	200 $\mu\text{g/mL}$ STD + H_2O	Vol (μL)	Glycogen ($\mu\text{g/mL}$)
1	200 μL + 0 μL	200	200
2	150 μL + 50 μL	200	150
3	100 μL + 100 μL	200	100
4	50 μL + 150 μL	200	50
5	0 μL + 200 μL	200	0

Transfer 10 μL standard and samples into separate wells of a clear flat-bottom microplate. If the sample contains glucose, transfer an additional 10 μL sample to another well for the Sample Blank.

3. **Working Reagent.** For each reaction well, mix 90 μL Assay Buffer, 1 μL Enzyme A, 1 μL Enzyme B and 1 μL Dye Reagent. For Sample Blank wells, prepare **Blank Working Reagent** by mixing 90 μL Assay Buffer, 1 μL Enzyme B, and 1 μL Dye Reagent (**No Enzyme A**). Transfer 90 μL Working Reagent into each Standard and Sample well. Transfer 90 μL Blank Working Reagent to each Sample Blank well. Tap plate to mix.

4. Incubate 30 min at room temperature. Read optical density at 570 nm (550-585 nm).

Fluorimetric Procedure:

For fluorimetric assays, the linear detection range is 0.2 to 20 $\mu\text{g/mL}$ glycogen. Follow steps 1-3 of the colorimetric procedure, but prepare 0, 5, 10, 15 and 20 $\mu\text{g/mL}$ Standard and use a black flat-bottom microplate. Incubate 30 min at room temperature and read fluorescence at $\lambda_{\text{ex}} = 530$ nm and $\lambda_{\text{em}} = 585$ nm.

CALCULATION

Subtract Blank reading ($\text{OD}_{570\text{nm}}$ or fluorescence intensity) from the standard reading values and plot the ΔOD or ΔF against standard concentrations. Determine the slope and calculate the glycogen concentration of the sample.

$$\text{Glycogen} = \frac{R_{\text{SAMPLE}} - R_{\text{BLANK}}}{\text{Slope}} \mu\text{g/mL}$$


R_{SAMPLE} and R_{BLANK} are the $\text{OD}_{570\text{nm}}$ or fluorescence intensity values of the sample and blank (water, or sample blank).

GENERAL CONSIDERATIONS

1. This assay is based on a kinetic reaction, the use of a multi-channel pipettor for adding the working reagent is recommended.
2. SH-group containing reagents (e.g., DTT, β -mercaptoethanol) may interfere with this assay and should be avoided in sample preparation.

MATERIALS REQUIRED, BUT NOT PROVIDED

Pipetting devices, centrifuge tubes, clear flat-bottom uncoated 96-well plates (e.g. VWR cat# 82050-760), optical density plate reader; black flat-bottom uncoated 96-well plates (e.g. VWR cat# 82050-676), fluorescence plate reader.



PUBLICATIONS

1. Hunter, A. L., et al (2020). Nuclear receptor REVERB α is a state-dependent regulator of liver energy metabolism. *Proceedings of the National Academy of Sciences*, 117(41), 25869-25879.
2. Primassin, S., et al (2011). Hepatic and muscular effects of different dietary fat content in VLCAD deficient mice. *Mol Genet Metab* 104(4): 546-551.
3. Yin, M et al (2011). Metformin improves cardiac function in a nondiabetic rat model of post-MI heart failure. *Am J Physiol Heart Circ Physiol* 301(2):H459-468.