

EnzyChrom™ Malate Dehydrogenase Assay Kit (EMDH-100)

Quantitative Colorimetric Kinetic Malate Dehydrogenase Activity Determination

DESCRIPTION

MALATE DEHYDROGENASE (MDH) (EC 1.1.1.37) is an enzyme which reversibly catalyzes the oxidation of L-malate to oxaloacetate in the presence of NAD. There are 2 isoforms in eukaryotic cells: MDH1 and MDH2. MDH1 found in the cytoplasm and plays a key part in the malate-aspartate shuttle for transporting malate into the mitochondria. MDH2 is a mitochondrial enzyme which participates in the TCA cycle that reversibly converts L-malate into oxaloacetate. Higher MDH activities are found in some neurodegenerative diseases such as Alzheimer's disease.

BioAssay Systems' non-radioactive, colorimetric MDH assay is based on the reduction of the tetrazolium salt MTT in a NADH-coupled enzymatic reaction to a reduced form of MTT which exhibits an absorption maximum at 565 nm. The increase in absorbance at 565 nm is proportional to the enzyme activity.

KEY FEATURES

Fast and sensitive. Linear detection range (20 µL sample): 0.5 to 65 U/L for 20 min reaction at 37°C.

Convenient and high-throughput. Homogeneous "mix-incubate-measure" type assay. Can be readily automated on HTS liquid handling systems for processing thousands of samples per day.

APPLICATIONS

MDH activity determination in biological samples (e.g. plasma, serum, erythrocytes, tissue and culture media.)

KIT CONTENTS (100 TESTS IN 96-WELL PLATES)

Assay Buffer: 10 mL	Enzyme A: 120 µL
NAD/MTT: 1 mL	Enzyme B: 120 µL
Substrate: 600 µL	Calibrator: 1.5 mL

Storage conditions. The kit is shipped on ice. Store all components at -20°C upon receiving. Shelf life: 6 months after receipt.

Precautions: reagents are for research use only. Normal precautions for laboratory reagents should be exercised while using the reagents. Please refer to Material Safety Data Sheet for detailed information.

PROCEDURES

This assay is based on a kinetic reaction. To ensure identical incubation time, addition of Working Reagent to samples should be quick and mixing should be brief but thorough. Use of a multi-channel pipettor is recommended. Assays can be executed at any desired temperature (e.g. 25°C or 37°C).

Sample Preparation: Serum and plasma are assayed directly.

Tissue: prior to dissection, rinse tissue in phosphate buffered saline (pH 7.4) to remove blood. Homogenize tissue (50 mg) in ~200 µL cold 50 mM potassium phosphate buffer, pH 7.5. Centrifuge at 14,000 × g for 10 min at 4°C. Remove supernatant for assay.

Cell Lysate: collect cells by centrifugation at 2,000 × g for 5 min at 4°C. For adherent cells, do not harvest cells using proteolytic enzymes; rather use a rubber policeman. Homogenize or sonicate cells in an appropriate volume of cold buffer containing 50 mM potassium phosphate (pH 7.5). Centrifuge at 14,000 × g for 10 min at 4°C. Remove supernatant for assay.

All samples can be stored at -20 to -80°C for at least one month.

Reagent Preparation: Equilibrate reagents to desired reaction temperature (37°C is recommended). Briefly centrifuge tubes before use.

Assay Procedure:

1. Transfer 100 µL H₂O (OD_{H2O}) and 100 µL Calibrator (OD_{CAL}) solution into wells of a clear flat bottom 96-well plate.
2. Transfer 20 µL H₂O into one well, this will be the blank. Transfer 20 µL of each sample into separate wells.

3. Prepare enough Working Reagent (WR) for all reaction wells by mixing, for each 96-well assay, 74 µL Assay Buffer, 8 µL NAD/MTT, 5 µL Substrate, 1 µL Enzyme A, 1 µL Enzyme B.

Add 80 µL WR to all samples and blank wells. Tap plate briefly to mix.

4. Read OD_{565nm} at time 10 min (OD₁₀) and time 30 min (OD₃₀) on a plate reader.

CALCULATION

Subtract the OD₁₀ from OD₃₀ for each sample to compute the ΔOD_S values, do the same for the blank to compute ΔOD_B. MDH activity can then be calculated as follows:

$$\text{MDH Activity} = \frac{\Delta\text{OD}_S - \Delta\text{OD}_B}{\epsilon_{\text{mtt}} \cdot l} \times \frac{\text{Reaction Vol } (\mu\text{L})}{t \text{ (min)} \cdot \text{Sample Vol } (\mu\text{L})}$$

$$= \frac{273}{t \text{ (min)}} \times \frac{\Delta\text{OD}_S - \Delta\text{OD}_B}{\text{OD}_{\text{CAL}} - \text{OD}_{\text{H}_2\text{O}}} \times n \text{ (U/L)}$$


where ϵ_{mtt} is the molar absorption coefficient of reduced MTT. l is the light path length which is calculated from the calibrator. OD_{CAL} and OD_{H2O} are OD_{565nm} (OD₁₀) values of the Calibrator and water. t is the difference in time between readings (20 min is the recommended time at 37°C). Reaction Vol and Sample Vol are 100 µL and 20 µL, respectively. n is the dilution factor if the sample needed to be diluted.

Unit definition: 1 Unit (U) of MDH will catalyze the conversion of 1 µmole of oxaloacetate and NADH per minute at pH 7.5.

Note: If sample MDH activity exceeds 65 U/L, dilute samples in water and repeat the assay. For samples with MDH activity < 1 U/L, the incubation time can be extended to 2 hours.

MATERIALS REQUIRED, BUT NOT PROVIDED

Pipetting devices and accessories (e.g. multi-channel pipettor), clear flat-bottom 96-well plates (e.g. VWR cat# 82050-760), centrifuge tubes and plate reader.



PUBLICATIONS

1. Venkat, S. et al (2017). Genetically encoding thioacetyl-lysine as a non-deacetylable analog of lysine acetylation in Escherichia coli. FEBS open bio, 7(11), 1805-1814.
2. Gan, Q et al (2016). Expanding the genetic code of Salmonella with non-canonical amino acids. Sci Rep 6, 39920.
3. Trivedi, AK et al (2015). Adaptation of oxidative phosphorylation to photoperiod-induced seasonal metabolic states in migratory songbirds. Comp. Biochem. Physiol. A 184: 34-40.

