

# Anti-Cytokeratin, Basic (Type II or HMW) (Epithelial Marker) Monoclonal Antibody

Catalog Number: M12092-1

#### **About KRT76**

This monoclonal antibody recognizes basic (Type II or HMW) cytokeratins, which include 67kDa (CK1); 64kDa (CK3); 59kDa (CK4); 58kDa (CK5); 56kDa (CK6); 52kDa (CK8). Twenty human keratins are resolved with two-dimensional gel electrophoresis into acidic (pI 6.0) subfamilies. The acidic keratins have molecular weights (MW) of 56.5, 55, 51, 50, 50', 48, 46, 45, and 40kDa. This monoclonal antibody recognizes the 65-67, 64, 59, 58, 56, and 52kDa keratins of basic subfamily. Many studies have shown the usefulness of keratins as markers in cancer research and tumor diagnosis. SPM115/SPM116 is a broad spectrum anti pan-keratin antibody cocktail, which differentiates epithelial tumors from non-epithelial tumors e.g. squamous vs. adenocarcinoma of the lung, liver carcinoma, breast cancer, and esophageal cancer.

#### **Overview**

Product Name	Anti-Cytokeratin, Basic (Type II or HMW) (Epithelial Marker) Monoclonal Antibody
Reactive Species	Chicken, Cow, Dog, Human, Monkey, Mouse, Rabbit, Rat
Description	Boster Bio Anti-Cytokeratin, Basic (Type II or HMW) (Epithelial Marker) Monoclonal Antibody (Catalog # M12092-1). Tested in Flow Cytometry, IF, WB, IHC applications. This antibody reacts with Human, Monkey, Cow, Dog, Rabbit, Mouse, Rat, Chicken.
Conjugate	Biotin
Application	Flow Cytometry, IF, IHC, WB
Clonality	Monoclonal Clone: SPM116
Formulation	Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.
Storage Instructions	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Host	Mouse
Uniprot ID	Q01546

#### **Technical Details**

Immunogen	Human epidermal keratin
Predicted Reactive Species	Pig, Rabbit
Cross Reactivity	Does not cross-react with primate, avian or amphibian GR.
Isotype	IgG1, kappa






888-466-3604 | support@bosterbio.com | www.bosterbio.com

Form	Liquid
Concentration	Purified antibody with BSA and azide at 200ug/ml
Purification	200ug/ml of antibody purified from Bioreactor Concentrate by Protein A/G.
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.  If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.  Some PubMed article(s) citing the expression level of this target are as follows:  Boster Bio's internal QC testing used: Flow Cytometry (1-2ug/million cells) Immunofluorescence (1-2ug/ml) Western Blot (1-2ug/ml f) Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 min at RT)(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95&degC followed by cooling at RT for 20 minutes) Optimal dilution for a specific application should be determined.



## Anti-Cytokeratin, Basic (Type II or HMW) (Epithelial Marker) Monoclonal Antibody (M12092-1) Images



Formalin-paraffin human Breast Carcinoma stained with Anti-Cytokeratin, HMW MAb (SPM116).

### Submit a product review to Biocompare.com





Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-Cytokeratin, Basic (Type II or HMW) (Epithelial Marker) Monoclonal Antibody