

pMX-STAT5A(1*6) Retroviral Vector

CATALOG NUMBER: RTV-331**STORAGE:** -80°C**QUANTITY AND CONCENTRATION:** 100 µL of bacterial glycerol stock

Background

Retroviruses are efficient tools for delivering heritable genes into the genome of dividing cells. Cell Biolabs' pMX retroviral vector is based on Moloney murine leukemia virus (MMLV). The viral *env* gene, produced by the package cell line, encodes the envelop protein, which determines the viral infectivity range. Transfection into a package cell line produces high-titer, replication-incompetent viruses. In addition to transfer and expression of exogenous genes in mammalian cells, recently, retroviruses have been used to express silencing RNAs (siRNA) to decrease the expression of target genes both *in vitro* and *in vivo*.

STAT5A and STAT5B are encoded by different genes; however, they share >90% amino acid identity. While initially referred to as mammary gland factor because of the role in mediating the effects of prolactin on the expression of β -casein in the mammary gland, mammary gland factor was renamed STAT5 when the cDNA was cloned and shown to be a member of the STAT family of transcription factors. Subsequently, two separate genes were identified and shown to be expressed in a wide range of tissues. STAT5A and STAT5B are now known to be differentially expressed in various tissues, and to have independent and distinct as well as common functions. STAT5B accounts for ~90% of the STAT5 in the liver; however, STAT5B and STAT5A are both required for the constitutive expression of certain GH-regulated liver cytochrome P450 enzymes.

The vector contains the ampicillin-resistance gene, MMLV LTRs, package signal and mouse constitutively active STAT5A(1*6) (H299R, S711F) (Figure 1).




Figure 1. Schematic representation of pMX-STAT5A(1*6) retroviral vector.

Safety Consideration

Remember that you will be working with samples containing infectious virus. Follow the recommended NIH guidelines for all materials containing BSL-2 organisms. Always wear gloves, use filtered tips and work under a biosafety hood.

References

1. Onishi, M., *et al.*, (1998) *Mol. Cell. Biol.* **18**, 3871-3879.
2. Morgenstern, J. P. and H Land. (1990) *Nuc. Acid Res.* **18**, 3587-3596.
3. Coffin, J. M. and H. E. Varmus, *Retroviruses*, Cold Spring Harbor Press, NY.
4. Onishi, M., Kinoshita, S., Morikawa, Y., *et al.*, (1996) *Exp. Hematol.* **24**, 324-327.
5. Kitamura T., *et al.*, (2003) *Exp. Hematol.* **31**, 1007-1014.

Recent Product Citation

1. Yu, Y. *et al.* (2012). Bcl11a is Essential for Lymphoid Development and Negatively Regulates p53. *J. Exp. Med.* **209**:2467-2483.

License Information

This product is licensed from Dr. Toshio Kitamura.

Warranty

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

This product is for RESEARCH USE ONLY; not for use in diagnostic procedures.

Contact Information

Cell Biolabs, Inc.
7758 Arjons Drive
San Diego, CA 92126
Worldwide: +1 858-271-6500
USA Toll-Free: 1-888-CBL-0505
E-mail: tech@cellbiolabs.com
www.cellbiolabs.com

©2012: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.