

Anti-Human Interferon Beta (MMHB-3) Neutralizing Monoclonal Antibody

Catalog Number: CM2015 **Size:** 50 µg

Description: Mouse Anti-Human Interferon Beta, neutralizing

Clone: MMHB-3

Concentration: 1.1 mg/ml

Buffer: phosphate buffered saline (PBS) containing 0.1% bovine serum albumin (BSA)

Antigen: human interferon beta

Isotype: IgG1 κ

Bioactivity: neutralizes human interferon beta; binds to human interferon beta with high affinity

Assay Used to Measure Bioactivity:

One neutralization unit is the amount of antiserum required to neutralize one unit of human interferon beta (hIFN-β) to a 50% endpoint. Interferon was titrated with the use of the cytopathic effect inhibition assay as described [Rubinstein, S., Familletti, P.C., and Pestka, S. (1981) "Convenient Assay for Interferons," *J. Virol.* **37**, 755-758; Familletti, P.C., Rubinstein, S., and Pestka, S. (1981) "A Convenient and Rapid Cytopathic Effect Inhibition Assay for Interferon," in *Methods in Enzymology*, Vol. 78 (S. Pestka, ed.), Academic Press, New York, 387-394]. In this antiviral assay for interferon about 1 unit/ml of interferon is the quantity necessary to produce a cytopathic effect of 50%. The units are determined with respect to the international reference standard for hIFN-β provided by the National Institutes of Health [see Pestka, S. (1986) "Interferon Standards and General Abbreviations," in *Methods in Enzymology* (S. Pestka, ed.), Academic Press, New York **119**, 14-23].

Usage: Western blot; immunoprecipitation; immunohistochemistry; ELISA

Shipping Temperature: Dry ice

Physical State of Product During Shipping: Frozen

Special Conditions/Comments: After receipt, this product should be kept at -20°C or below for retention of full activity. When thawing, the contents of the tube should be apportioned in separate tubes so that freezing and thawing is kept to a minimum.

Usage: **For Research Use Only. Not for use in diagnostic procedures.**

Country of Origin: USA

* This is a generic data sheet and may be subject to change. Please see the package insert shipped with your product for current data.



Cell Sciences, Inc.
480 Neponset Street
Bldg 12A
Canton, MA 02021

Toll Free: 888-769-1246
Phone: 781-828-0610
Fax: 781-828-0542

E-mail: info@cellsciences.com
Web Site: www.cellsciences.com