

Anti-Human Interferon Beta

Polyclonal Antibody

Catalog Number: CP2006 Size: 1 vial of 2 x 10⁴ units in 0.215 ml

Lot Number: 1881

Description: Sheep Anti-Human Interferon Beta

Bioactivity: $9.3 \times 10^4 \text{ units/ml}$ **Concentration:** $9.3 \times 10^4 \text{ units/ml}$

Form: neat serum

Antigen: recombinant human interferon beta (grown in *E. coli*)

Assay Used to Measure Bioactivity:

One neutralization unit is the amount of antiserum required to neutralize one unit of human interferon beta (hIFN- β) to a 50% endpoint. Interferon was titrated with the use of the cytopathic effect inhibition assay as described [Rubinstein, S., Familletti, P.C., and Pestka, S. (1981) "Convenient Assay for Interferons," *J. Virol.* **37**, 755-758; Familletti, P.C., Rubinstein, S., and Pestka, S. (1981) "A Convenient and Rapid Cytopathic Effect Inhibition Assay for Interferon," in *Methods in Enzymology*, Vol. 78 (S. Pestka, ed.), Academic Press, New York, 387-394]. In this antiviral assay for interferon about 1 unit/ml of interferon is the quantity necessary to produce a cytopathic effect of 50%. The units are determined with respect to the international reference standard for hIFN- β provided by the National Institutes of Health [see Pestka, S. (1986) "Interferon Standards and General Abbreviations," in *Methods in Enzymology* (S. Pestka, ed.), Academic Press, New York **119**, 14-23]. This material is prepared specifically for effective neutralization of hIFN- β .

Usage: Western blot; immunoprecipitation; immunohistochemistry; ELISA

Special Conditions/Comments: After receipt, this product should be kept at -20°C or below for

retention of full activity. When thawing, the contents of the tube should be apportioned in separate tubes so that freezing and

E-mail: info@cellsciences.com

Web Site: www.cellsciences.com

thawing is kept to a minimum.

Toll Free: 888-769-1246

Phone: 781-828-0610

Fax: 781-828-0542

Usage: For Research Use Only. Not for use in diagnostic procedures.

Country of Origin: USA