

Human Chorionic Gonadotropin, ß-HCG, clone 2092 Monoclonal Antibody

Catalog No.: MON2001-1

Quantity: 1 ml

Specificity

The antibody reacts with high specificity with ß-HCG.

Cross reactivity with: FSH: not detectable. LH: slight detectable. MSH: not detectable.

Immunoglobulin type

Murine IgG₁

Use

The antibody stains HCG containing cells in frozen or paraffin embedded tissue sections with immunohistochemistry and immunofluorescence methods and reacts with HCG in ELISA tests. It cannot be used for RIA with radio-labeled antigen.

Instructions for use

For use in staining of tissue sections the antibody has to be diluted 1:5-1:25, preferably in phosphate buffered saline.

Presentation

Each vial contains 1 ml tissue culture supernatant with approximately 20 μ g/ml antibody and 1% BSA and 10mM sodium azide in PBS. Sufficient for 50-250 tests. Cat. nr. MON2001-1.

Each vial contains 5 ml tissue culture supernatant with approximately 20 μ g/ml antibody and 1% BSA and 10mM Sodium azide in PBS. Sufficient for 250-1250 tests. Cat.nr. MON2001-5.

Toll Free: 888-769-1246

Phone: 781-828-0610

Fax: 781-828-0542

E-mail: info@cellsciences.com

Web Site: www.cellsciences.com

FOR RESEARCH USE ONLY, NOT FOR DRUG, DIAGNOSTIC OR OTHER USE.