

CD106, VCAM-1, Biotin, clone 1G11B1 Monoclonal Antibody

Catalog No.: MON 6017B

Quantity: 0.5 ml

Specificity

The antibody reacts with VCAM-1, the vascular cell adhesion molecule (1), which is a member of the immunoglobulin superfamily, C2 subset. The antigen is induced on human endothelial cells by $\text{TNF}\alpha$, IL-1, IFN- γ or endotoxin. The ligands for VCAM-1 are integrin α_4 , β_1 (VLA-4) and integrin α_4 , β_7 , integrins that are constitutively expressed by thymocytes, lymphocytes and monocytes.

Immunoglobulin type

IgG₁

Use

The antibody is extremely useful for staining of VCAM-1 expressing endothelial cells. It permits staining of in vitro cultured cells and frozen tissue sections. It inhibits cellular adhesion. this antibody is useful in ELISA techniques.

Instructions for use

In vitro cultured cells can be fixed with 1% paraformaldehyde. Tissue sections are advised to be fixed for 10 min. in pure acetone and followed by incubation for 10 min in chloroform. Incubation with a pretested dilution of the antibody is advised to be incubated with an enzyme (alkaline phosphatase) conjugated streptavidin. For selection of the most useful dilution in a given situation a test staining with cells or tissue known to express the antigen should be performed. To this end either cultured endothelial cells or a small fresh skin biopsy can be incubated for 6 hours with $\text{TNF}\alpha$ (1 ng/ml), IL-1 (100 u/ml) or LPS (1 $\mu\text{g/ml}$) in tissue culture medium at 37°C. As negative control it is advised to use as first antibody a control murine IgG₁ antibody.

Presentation

0.5 ml. containing approximately 100 $\mu\text{g/ml}$ purified Ig with 0.02% sodium azide, sufficient for staining of 100 frozen tissue samples.

Literature:

1. Polte T., et al., 1990, N.A.R. 15, 5901.
2. Thornhill, M., et al., 1991, J. Immunol. 146, 592.
3. Kyan-Aung, U., et al., 1991, Am. Review. 5, 445-450.

FOR RESEARCH USE ONLY, NOT FOR DRUG, DIAGNOSTIC OR OTHER USE.



Cell Sciences, Inc.
480 Neponset Street
Bldg 12A
Canton, MA 02021

Toll Free: 888-769-1246
Phone: 781-828-0610
Fax: 781-828-0542

E-mail: info@cellsciences.com
Web Site: www.cellsciences.com