

## Rabbit Anti-Cat NTS/Neurotensin pAb

<b>Catalog No.</b>	PS224	<b>Quantity:</b>	50 µl
<b>Specificity:</b>	Positive control: Bouin-fixed paraffin or Stefanini-fixed frozen sections of cat ileum. Absorption with 10-100 µg neurotensin per ml diluted antiserum abolishes the staining while glucagon, insulin, gastrin, GIP and VIP do not.		
<b>Host:</b>	Rabbit		
<b>Immunogen:</b>	Synthetic bovine neurotensin, conjugated to BSA. Tested on cat tissues.		
<b>Formulation:</b>	50 µl undiluted, freeze-dried antiserum. Dissolve the antiserum in 50-100 µL distilled water, and dilute further in 0.1M PBS with 1% BSA and 0.1% NaN <sub>3</sub> .		
<b>Applications:</b>	Works on: Paraffin and frozen sections. 1:2000-1:3000 using PAP, or 1:1000-1:2000 using FITC with overnight incubation at +4°C.		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

