

Mouse TNF-alpha ELISPOT Kit with Sterile PVDF Plates

Catalog No: CKM032A
CKM032B

Lot: xxxxxxxx


Size: 2-plate format
5-plate format

Expiration Date: xx.xx.xxxx

Contents of the kit:

Catalog #	Item	Storage Temp.	Quantity (2-plate format)	Quantity (5-plate format)
CKM032-A	Coating antibody (lyophilized)	4°C	1 vial	1 vial
CKM032-B	Biotinylated detection antibody (lyophilized)	4°C	1 vial	1 vial
CKM032-C	Streptavidin-HRP conjugate (lyophilized)	≤-20°C	1 vial	1 vial
	AEC coloring system:			
CKM032-D	I. AEC stock solution	≤-20°C (in dark)	1 vial (4 ml)	1 vial (4 ml)
CKM032-E	II. Substrate buffer capsule	RT	1 vial (2.5 ml)	1 vial (5 ml)
CKM032-F	Blocking stock solution (10x)	4°C	1 vial (4 ml)	1 vial (10 ml)
CKM032-G	Dilution buffer R (10x)	4°C	1 vial (4 ml)	1 vial (10 ml)
CKM032-H	Tween-20	RT	1 vial (5 ml)	1 vial (5 ml)
CKM032-P	Sterile 96-well PVDF plate (Millipore MSIPS4510)	RT	2	5
CKM032-S	Adhesive cover slip	RT	5	10

Typical data for CKM032



Stimulus: none (1×10^5 cells/well)



Stimulus: Concanavalin A (1×10^5 cells/well)

Example of TNF- α specific spots produced by mouse (CD-1) splenocytes.
Preincubation: 24h, incubation ELISPOT: 20h.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences, Inc.
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Web Site: www.cellsciences.com

cellsciences.com



Cell Sciences, Inc.
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Web Site: www.cellsciences.com